

THERMO JACKET SET

Product Name SV-120
Category CATEGORY III
N° Certificate 63-20-07-R01
Date of Certificate 05.04.2021-05.12.2025

PRODUCT INFORMATION

Color : Blue, White, Green
Packaging : Single pack of 60 pcs.
Sizes : S-4XL
Material : 20-80 GR
Classification : SS-SMS
Product Type : THERMO JACKET SET
Model Definition : THERMO JACKET SET



DESCRIPTION SMS SURGICAL GOWN SET

Surgical suits are designed for both men and women. Surgical suits can be easily used by doctors, nurses, anesthesia technicians, midwives and other health personnel. It is 2 different products and 1 set consisting of upper and lower parts. The lower and upper parts are produced from the same color and laminated fabric. It is made of sweat-proof and stain-resistant fabric. Our surgical jacket sets have buttons that can be opened easily. In this way, it provides more comfortable use in putting on and taking off.

PPE Specification

Classification
Abrasion Resistance
Tear Resistance
Tensile Strength
Puncture Resistance
Liquid Repellency
Resistance to Liquid Penetration
Seam Strength
Wet Bacterial Penetration
Dry Microbial Penetration
Phi-x174 Bacteriophage

Performance Levels

: Type PB 3-B, Type PB 4-B, Type PB 6-B
: 1
: 2
: 1
: 1
: NaOH:3, H₂SO₄: 3
: NaOH:3, H₂SO₄: 3
: 2
: 1
: 1
: 3

THERMO JACKET SET



M E D İ K A L

Product Name SV-120
Category CATEGORY III
N° Certificate 63-20-07-R01
Date of Certificate 05.04.2021-05.12.2025

IDENTIFICATION

EN ISO 13688 Protective clothing – General requirements
EN 13034 + A1 Protective clothing Against Liquid Chemicals (Type PB-6-B)
EN 14126 Protective clothing Performance requirements and tests methods for protective clothing against infective agents
ISO 9001: 2015 Quality Management System
ISO 10002: 2018 Customer Satisfaction Management System
ISO 13485: 2016 Medical Devices Management System
ISO 14001:2015 Environment Management System
ISO 45000:2018 Occupational and Safety Management System
FDA

WARNINGS

This is an indicative guide only. It must not be used as the only method for selecting the protective garment. Before using any protective clothing. The user must have read and understood all the instructions for use of the product. It is necessary to observe the legislation in force in each individual country. The selection of the most appropriate PPE or accessory depends on the particular situation and should only be made by a competent person who is aware of the specific working conditions and the limitations of PPE and accessories. It is the employer's responsibility to define the suitability of these products for the intended use. This information is subject to review at any time.



M E D İ K A L

PODIMA MEDİKAL VE TEKSTİL SAN.TİC. LTD. STİ.
Y.DUDULLU MAH. BAYRAK CAD. NO:30 D:130
UMRANIYE / İSTANBUL 34775

THERMO JACKET SET

PODiMA

M E D İ K A L

Product Name SV-120
Category CATEGORY III
N° Certificate 63-20-07-R01
Date of Certificate 05.04.2021-05.12.2025

PACKAGING SIZE



Width*Depth*Height

- 1) Carton Size **60*40*40 cm**
- 2) Total Pallet Size **120*80*255 cm** (24 cartons in 1 pallet)

1 Container : 852 CTN (Without pallets)

1 Container : 768 CTN in 32 Pallets

- Single pack of 60 pcs.



PODiMA

M E D İ K A L

PODiMA MEDİKAL VE TEKSTİL SAN.TİC. LTD. STİ.
Y.DUDULLU MAH. BAYRAK CAD. NO:30 D:130
UMRANIYE / İSTANBUL 34775