

PRODUCT INFORMATION

| | |
|----------------|---------------------------------------|
| COLORS | : Blue, White, Yellow, Green |
| PACKAGING | : From 1 pcs/bag to 10 pcs/bag |
| SIZES | : S-4XL |
| CUFFS | : Knitted / Elastic |
| NECK CLOSURE | : Velcro Fastener / Ties |
| BELT CLOSURE | : Standard / Ties |
| SEWING TYPE | : Overlock / Heat Sealed / Ultrasonic |
| STERILIZATION | : Non-Sterile / Sterile (EO) |
| MATERIAL | : SMS |
| WEIGHT | : From 20 GSM to 80 GSM |
| CLASSIFICATION | : CAT III |



Please contact us for products technical file. Using or copying without our permission will have legal consequences.

Ürün teknik dosyası için lütfen bizimle iletişime geçiniz. İznimiz olmadan kullanılması veya kopyalanması hukuki sonuçlar doğuracaktır.