



High protection
made with
comfort

Protective Clothing

Your reliable partner for personal protective equipment, including coveralls, protective gowns, hoods, shoe covers, nitrile gloves and protective masks.



STANDARD
100



Gowns

Disposable visitor, examination and surgical gowns for medical sector to protect health professionals, patients and visitors against viruses and germs. Also applicable in food industry, manufacturing, laboratories and pharmacies.

Coveralls

Disposable coveralls ensure protection against infective agents, chemical splashes and other particles. Manufactured of a soft and breathable material coveralls offer comfort and full-body protection.



Hoods

Disposable hoods ensure multiple protections for open wounds during surgeries. Additionally, disposable hoods avoid food to be contaminated by hair or dandruff in food industry or gastronomy.

Shoe covers

Applicable in medical sector and food industry, shoe covers avoid contamination and offer protection both for ground (e.g. sandy shoes on damp ground) and shoes (e.g. protection from dust or paints during renovation works).





Towel & Sheets



CONTENTS

Gowns	4
Coveralls	19
Hoods & Bonnets	26
Shoe Covers	28
Mask	31
Bed Sheet	34
Towels & Sheets	35
Certificates	41
Production Process	42

CERTIFICATIONS AND STANDARDS



Gowns

FOR SAFETY OF ALL IN HEALTHCARE SETTINGS

Examination gowns are especially suitable for general medical practice and sensitive areas such as operating theatres. This gown is used to protect health professionals against contamination with bacteria and viruses in hospitals, operating theatres, laboratories or other medical facilities.

 **PODiMA**

 **SUVICOM®**

 **rowimed®**

 **Raddax®**



DISPOSABLE VISITOR GOWN

Visitor gowns are used for minimal risk situations which contain slight barrier to small fluid penetration. The gowns are breathable and comfortable to wear.

COLORS



SERVICES




SIZE	XS	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	56	56	56	56	56	56	58	58	+/- 3 cm
B	58	59	60	61	62	63	64	65	+/- 3 cm
C	105	110	115	120	125	130	135	140	+/- 3 cm
D	130	135	140	145	150	155	160	165	+/- 3 cm

SV- 180

DETAILS

Model: SV-180 series
Fabric: PP (spunbond+spunbond)
Weight: 17 - 80 GSM
Cuff: Elastic/Knitted
Color: White/Yellow/Blue/Green
Stitch: Overlock
Other: Waist and neck straps for a secure fit
Optional: Tie on neck - Velcro neck
Sterile/Non Sterile

For more informatin
please refer to the data sheets.



DISPOSABLE ISOLATION GOWN (PP+PE)

Examination gowns are especially suitable for general medical practice and sensitive areas such as operating theatres. This gown is used to protect health professionals against contamination with bacteria and viruses in hospitals, operating theatres, laboratories or other medical facilities.

COLORS



SERVICES



SIZE	XS	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	56	56	56	56	56	56	58	58	+/- 3 cm
B	58	59	60	61	62	63	64	65	+/- 3 cm
C	105	110	115	120	125	130	135	140	+/- 3 cm
D	130	135	140	145	150	155	160	165	+/- 3 cm

SV- 140

DETAILS

Model: SV- 140 series
 Fabric: PP + PE (polyethylene) laminated
 Weight: 30-80 GSM
 Cuff: Elastic/Knitted
 Color: Blue/Yellow/Green/White
 Stitch: Overlock / Ultrasonic seam
 Other: Waist and neck straps for a secure fit
 Tie on neck - Velcro neck
 Sterile/Non Sterile

For more information
 please refer to the data sheets.



Velcro neck



Tie on belt



Knitted cuff



SSMMS/ Comfort



DISPOSABLE SURGICAL GOWN (SMMS)

Surgical gowns are made of comfortable non-woven fabric and contain high-quality four layers of SMMS (spunbond+meltblown + meltblown + spunbond) fabric. This gown is used to protect health professionals against contamination with bacteria and viruses in hospitals, operating theatres, laboratories or other medical facilities.

SV- 141

DETAILS

Model: SV- 141 series
Fabric: SMMS (spunbond+meltblown +meltblown+spunbond) + Polypropylene non-woven (PP)
Weight: 25-80 GSM
Cuff: Elastic/Knitted
Color: White/Yellow/Blue/Green
Stitch: Overlock/Ultrasonic
Other: Waist and neck straps for a secure fit
Optional: Sterile/Non Sterile Velcro neck

For more informatin please refer to the data sheets.

DISPOSABLE SURGICAL GOWN (SSMMS)

Surgical gowns are made of comfortable non-woven fabric and contain high-quality four layers of SSMMS (spunbond +spunbond + meltblown + meltblown+spunbond) fabric. This gown is used to protect health professionals against contamination with bacteria and viruses in hospitals, operating theatres, laboratories or other medical facilities.

SV- 142

DETAILS

Model: SV- 142 series
Fabric: SSMMS (spunbond+spunbond+ meltblown+meltblown+spunbond + Polypropylene non-woven (PP)
Weight: 25-80 GSM
Cuff: Elastic/Knitted
Color: White/Yellow/Blue/Green
Stitch: Overlock/Ultrasonic
Other: Waist and neck straps for a secure fit
Optional: Sterile/Non Sterile Velcro neck


For more informatin please refer to the data sheets.

COLORS




SERVICES





Triplex Reinforced Breathable Laminated Fabric for Surgical Gown

- 1 SSMMS Spunbond
- 2 Hotmelt
- 3 Triple - Layered PE Film
- 4 Hotmelt
- 5 SSMMS Spunbond



HIGH PROTECTION SURGICAL GOWN (SSMLMS + PE)

Protective surgical gown is made of light and breathable material to use in healthcare environments to protect health professionals and staff against the transmission of microorganisms and body fluids. It is also comfortable and easy to tie at the waist and collar at the back.

COLORS



SERVICES



SIZE	XS	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	56	56	56	56	56	56	58	58	+/- 3 cm
B	58	59	60	61	62	63	64	65	+/- 3 cm
C	105	110	115	120	125	130	135	140	+/- 3 cm
D	130	135	140	145	150	155	160	165	+/- 3 cm

SV- 143

DETAILS

- Model: SV- 143 series
Fabric: 17 gsm SSMMS Nonwoven + 1,5 gsm Hotmelt + 18 gsm Triple-layered Breathable Medical PE Film + 1,5 gsm Hotmelt + 17 gsm SSMMS Nonwoven
Weight: 55 - 90 GSM
Cuff: Elastic/Knitted
Color: Blue
Stitch: Overlock/ Ultrasonic seam
Other: Waist and neck straps for a secure fit
Optional: Static/Antistatic
Tie on neck - Velcro neck
Sterile/Non Sterile

For more informatin
please refer to the data sheets.



DISPOSABLE Reinforced SURGICAL GOWN

Highly protective (Level 4) surgical gowns are used in healthcare environments, particularly operating theatres, to protect health professionals and staff against the transmission of microorganisms and body fluids. These gowns ensure highest resistance to moisture and bacteria penetration, puncture, tear, abrasion and toxic substances.

COLORS



SERVICES




SIZE	XS	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	56	56	56	56	56	56	58	58	+/- 3 cm
B	58	59	60	61	62	63	64	65	+/- 3 cm
C	105	110	115	120	125	130	135	140	+/- 3 cm
D	130	135	140	145	150	155	160	165	+/- 3 cm

SV- 144

DETAILS

- Model: SV- 144 series
- Fabric: SMS + PE (spunbond+meltblown +spunbond + polyethylene laminated)
- Weight: 45-80 GSM
- Cuff: Elastic/Knitted
- Color: Blue
- Stitch: Overlock/ Ultrasonic seam
- Other: Waist and neck straps for a secure fit
- Optional: Tie on neck - Velcro neck
- Sterile/Non Sterile

For more informatin
please refer to the data sheets.



DISPOSABLE **Reinforced** ISOLATION GOWN (PP+PE)

Examination gowns are especially suitable for sensitive areas such as operating theatres. It is particularly characterized by its tightness, therefore it is also certified to EN14126.

COLORS



SERVICES



SIZE	XS	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	56	56	56	56	56	56	58	58	+/- 3 cm
B	58	59	60	61	62	63	64	65	+/- 3 cm
C	105	110	115	120	125	130	135	140	+/- 3 cm
D	130	135	140	145	150	155	160	165	+/- 3 cm

SV- 144

DETAILS

Model: SV- 144 series
Fabric: PP + PE (polyethylene)
laminated (**partially reinforced**)
Weight: 35-80 GSM
Cuff: Elastic/Knitted
Color: Blue/Yellow/Green/White
Stitch: Overlock / Ultrasonic seam
Other: Waist and neck straps for a secure fit
Optional: Tie on neck - Velcro neck
Sterile/Non Sterile

For more information
please refer to the data sheets.



DISPOSABLE Reinforced SURGICAL GOWN / Ultimate

Highly fluid repellent gown for advanced protection. Ultimate offers maximum breathability for ultimate comfort.
High performance gown with impermeable front and sleeves for added protection against strike through. The gown cuff is compatible with gloves, preventing gloves from sliding off.
Available in High and Standard performance gowns. Conformable material. Shoulder pleats allow freedom of movement for the upper body.

COLORS



SERVICES



SV- 144

DETAILS

- Model: SV- 144 series
Fabric: SMS + PE (spunbond+meltblown +spunbond + polyethylene laminated)
Weight: 45-80 GSM
Cuff: Elastic/Knitted
Color: Blue
Stitch: Overlock/ Ultrasonic seam
Other: Waist and neck straps for a secure fit
Optional: Tie on neck - Velcro neck
Sterile/Non Sterile



SIZE	XS	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	56	56	56	56	56	56	58	58	+/- 3 cm
B	58	59	60	61	62	63	64	65	+/- 3 cm
C	105	110	115	120	125	130	135	140	+/- 3 cm
D	130	135	140	145	150	155	160	165	+/- 3 cm

For more information,
please refer to the data sheets.



DISPOSABLE EXAMINATION GOWN (SMS)

Disposable Surgical gowns are made of soft “spunbond + meltblown + spunbond” fabric. The gown is robust, breathable and liquid resistant with bacteria barrier and antiallergic protection. The gown also ensures protection against blood and other body fluids.

COLORS



SERVICES



SIZE	XS	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	56	56	56	56	56	56	58	58	+/- 3 cm
B	58	59	60	61	62	63	64	65	+/- 3 cm
C	105	110	115	120	125	130	135	140	+/- 3 cm
D	130	135	140	145	150	155	160	165	+/- 3 cm

SV- 145

DETAILS

- Model: SV- 145 series
- Fabric: SMS (spunbond+meltblown+spunbond)
- Weight: 20-80 GSM
- Cuff: Elastic/Knitted
- Color: White/Yellow/Blue/Green
- Stitch: Overlock - Ultrasonic
- Other: Waist and neck straps for a secure fit
- Optional: Tie on neck - Velcro neck
- Sterile/Non Sterile

For more informatin
please refer to the data sheets.




DISPOSABLE SURGICAL GOWN
(Thermoband Striped Gown)

Protective gown is made of light and breathable material to use in clinical and laboratory environments to protect health professionals, staff and visitors against the transmission of microorganisms and body fluids.

COLORS



SERVICES



SIZE	XS	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	56	56	56	56	56	56	58	58	+/- 3 cm
B	58	59	60	61	62	63	64	65	+/- 3 cm
C	105	110	115	120	125	130	135	140	+/- 3 cm
D	130	135	140	145	150	155	160	165	+/- 3 cm

SV- 146

DETAILS

- Model: SV- 146 series
Fabric: PP + PE +Thermoband
Weight: 40 - 80 GSM
medical PE film
Cuff: Elastic/Knitted
Color: Blue/Yellow/Green/White
Stitch: Overlock/Welded
Type: PB-3B / PB-4B / PB- 6B AAMI Level 2 - 3
Other: Waist and neck straps for a secure fit
Optional: Sterile/Non Sterile
Velcro neck
Ultrasonic seam

For more informatin
please refer to the data sheets.



Raglan Sewing



Standard Sewing

DISPOSABLE SURGICAL GOWN / Classic

Soft, breathable comfort and ease of movement in the process. Surgical gowns Classic are designed to protect you during most procedures.

Barrier surgical gown Classic is made from a soft and breathable spunlace material for extra comfort.

Specially treated to repel low-tension surface fluids.

Forward facing long sleeves are designed to fit your natural working position and reduce the risk of glove sliding.

COLORS



SERVICES



SIZE	XS	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	56	56	56	56	56	56	58	58	+/- 3 cm
B	58	59	60	61	62	63	64	65	+/- 3 cm
C	105	110	115	120	125	130	135	140	+/- 3 cm
D	130	135	140	145	150	155	160	165	+/- 3 cm

SV- 147

DETAILS

Model:	SV- 147 series
Fabric:	Spunlace / Sontara
Weight:	40 - 95 GSM
Cuff:	Elastic/Knitted
Color:	Blue/White
Stitch:	Overlock/Welded
Type:	PB-3B / PB-4B / PB- 6B AAMI Level 2 - 3
Other:	Waist and neck straps for a secure fit
Optional:	Sterile/Non Sterile
	Velcro neck
	Ultrasonic seam

For more information,
please refer to the data sheets.




V-NECK SET

V-neck set offer long-lasting warmth for the staff of operating room. The jacket is equipped with pockets based on requests


COLORS



SERVICES



	S	M	L	XL	2XL	3XL	Avarage
A	74	75	76	77	78	79	+/- 3cm
B	60	62	64	66	68	70	+/- 3cm
C	35	36	37	38	39	40	+/- 3cm
D	18	19	20	21	22	23	+/- 3cm



	S	M	L	XL	2XL	3XL	Avarage
A	102	104	106	108	110	112	+/- 3cm
B	52	54	57	61	64	67	+/- 3cm
C	54	56	59	63	66	69	+/- 3cm
D	21	22	23	24	25	26	+/- 3cm

SV- 120

DETAILS


Model: SV- 120 series
Fabric: SS-SMS-SMMS
Weight: 30-80 GSM
Cuff: Knitted
Color: Blue/Yellow/Green/White
Stitch: Overlock/Welded
Other: Desings with two or three pockets
Optional: Sterile/Non Sterile

For more informatin
please refer to the data sheets.



THERMAL SET & THERMAL JACKET

Thermal jackets offer long-lasting warmth for the staff of operating room. The jacket is equipped with pockets based on requests.



	S	M	L	XL	2XL	3XL	Avarage
A	73	75	77	79	81	83	+/- 3cm
B	55	57	59	61	63	65	+/- 3cm
C	59	60	61	62	63	64	+/- 3cm
D	28	28	28	28	29	29	+/- 3cm

	S	M	L	XL	2XL	3XL	Avarage
A	102	104	106	108	110	112	+/- 3cm
B	52	54	57	61	64	67	+/- 3cm
C	56	56	59	63	66	69	+/- 3cm
D	21	22	23	24	25	26	+/- 3cm

SV- 121

DETAILS

Model: SV- 121 series
Fabric: SS-SMS-SMMS
Weight: 30-80 GSM
Cuff: Knitted
Color: Blue/Yellow/Green/White
Stitch: Overlock/Welded
Other: Desings with two or three pockets
Optional: Sterile/Non Sterile

For more informatin
please refer to the data sheets.

SV- 122

DETAILS

Model: SV- 122 series
Fabric: SS-SMS-SMMS
Weight: 30-80 GSM
Cuff: Knitted
Color: Blue/Yellow/Green/White
Stitch: Overlock/Welded
Other: Desings with two or three pockets
Optional: Sterile/Non Sterile

For more informatin
please refer to the data sheets.

COLORS



SERVICES






KIMONO SET

Kimono Set are made of soft “spunbond + meltblown + spunbond” fabric. The gown is robust, breathable and liquid resistant with bacteria barrier and antiallergic protection. The gown also ensures protection against blood and other body fluids.


COLORS



SERVICES



	S	M	L	XL	2XL	3XL	4XL	Avarage
A	78	81	84	87	90	93	96	+/- 3cm
B	125	130	135	140	145	150	155	+/- 3cm
C	26	27	28	29	30	31	32	+/- 3cm





	S	M	L	XL	2XL	3XL	Avarage
A	102	104	106	108	110	112	+/- 3cm
B	52	54	57	61	64	67	+/- 3cm
C	54	56	59	63	66	69	+/- 3cm
D	21	22	23	24	25	26	+/- 3cm

SV- 123

DETAILS

- Model: SV- 123 series
- Fabric: SMS (spunbond+meltblown+spunbond)
- Weight: 20-80 GSM
- Cuff: Elastic/Knitted
- Color: White/Yellow/Blue/Green
- Stitch: Overlock - Ultrasonic
- Other: Waist and neck straps for a secure fit
- Optional: Tie on neck - Velcro neck
- Sterile/Non Sterile

For more informatin
please refer to the data sheets.



LAB COAT

Disposable Non - woven Lab Coat is made of non woven PP or SMS to seal the surface and provide protection from fluids. Non-toxic, no odor, does not irritate the skin .
The closure with velcro, zipper or buttons make the lab coat easy to take on or take off.
This breathable lab coat which is fit for use in food processing, dental clinics and science laboratories, and offers unmatched comfort and protection to the user's clothes from dry stuff.

COLORS



SERVICES



SIZE	S	M	L	XL	2XL	3XL	4XL	Avarage
A	96	99	102	105	108	111	114	+/- 3cm
B	61	64	67	70	73	76	79	+/- 3cm
C	56	56	56	56	56	58	58	+/- 3cm

SV- 189

DETAILS

Model: SV- 189 series
Fabric: SS-SMS-SMMS
Weight: 20-80 GSM
Cuff: Elastic, Knitted
Color: Blue/Yellow/Green/White
Stitch: Overlock/Welded
Other: Desings with two or three pockets
Optional: Sterile/Non Sterile

For more informatin
please refer to the data sheets.

Coveralls

FOR SAFETY OF ALL IN HEALTHCARE SETTINGS

Disposable coveralls ensure protection against infective agents, chemical splashes and other particles. Manufactured of a soft and breathable material coveralls offer comfort and full-body protection.

 **PODiMA**

 **SUVICOM®**

 **rowimed®**

 **Razdax®**




DISPOSABLE COVERALL TYPE 5B/6B

Disposable coverall ensures protection against infective agents, chemical splashes and other particles. Manufactured of a soft and breathable material this coverall offers comfort and high-level protection. The coverall is equipped with hood, self-adhesive zipper, secure and convenient closures with bands on hoods, arms and legs. These coveralls can be applicable in various areas such as handling chemicals in dust and powder form, varnishing and protection against paint splashes, fiberglass production and processing, industrial cleaning and maintenance, automotive industry, grinding and polishing, mining and quarrying, demolition and refurbishment, working with asbestos, wood and metal processing, construction industry and pharmaceutical industry.

COLORS



SERVICES



SIZE	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	55	60	65	70	75	80	85	+/- 3 cm
B	188	190	192	194	196	198	200	+/- 3 cm
C	74	76	78	80	82	84	86	+/- 3 cm

SV- 100

DETAILS

- Model: SV- 100 series
- Fabric: Spunbond hotmelt polyethylen, PP + PE
- Weight: 40-80 GSM
- Cuff: Elastic
- Color: White
- Stitch: Overlock
- Other: Waist and neck straps for a secure fit
- Optional: Static/Antistatic
Sterile/Non Sterile
Ultrasonic seam

For more informatin
please refer to the data sheets.




DISPOSABLE COVERALL (SMS - PP+PE) Without Hood

These coveralls offer comfort and are durable, breathable and water-resistant. The fabric is also silicone- and allergy-free. The coverall is equipped with hood, self-adhesive zipper, secure and convenient closures with bands on hoods, arms and legs. The welded stripes offer greater protection for liquids. These coveralls can be applicable in medical and industrial settings and conditions in which one works with dangerous biological and chemical substances.

COLORS



SERVICES



SIZE	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	55	60	65	70	75	80	85	+/- 3 cm
B	160	162	164	166	168	170	172	+/- 3 cm
C	74	76	78	80	82	84	86	+/- 3 cm

SV- 107

DETAILS

Model: SV- 107 series
Fabric: SMS, PP+PE
Weight: 20-80 GSM
Cuff: Elastic
Color: White
Stitch: Overlock
Optional: Static/Antistatic

For more informatin
please refer to the data sheets.




DISPOSABLE COVERALL TYPE 3B/4B THERMO BAND STRIPED

Disposable coverall ensures protection against infective agents, chemical splashes and other particles. Manufactured of a soft and breathable material this coverall offers comfort and protection. The coverall is equipped with hood, self-adhesive zipper, secure and convenient closures with bands on hoods, arms and legs. These coveralls can be applicable in various areas such as handling chemicals in dust and powder form, varnishing and protection against paint splashes, fiberglass production and processing, industrial cleaning and maintenance, automotive industry, grinding and polishing, mining and quarrying, demolition and refurbishment, working with asbestos, wood and metal processing, construction industry and pharmaceutical industry.

COLORS



SERVICES



SIZE	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	55	60	65	70	75	80	85	+/- 3 cm
B	188	190	192	194	196	198	200	+/- 3 cm
C	74	76	78	80	82	84	86	+/- 3 cm

SV- 110

DETAILS

- Model: SV- 110 series
- Fabric: Spunbond hotmelt polyethylen, PP + PE
- Weight: 55-80 GSM
- Cuff: Elastic
- Color: White
- Stitch: Ultrasonic stripe
- Other: Waist and neck straps for a secure fit
- Optional: Static/Antistatic
Sterile/Non Sterile

For more informatin
please refer to the data sheets.




DISPOSABLE COVERALL TYPE 3B/4B ULTRASONIC

Disposable coverall ensures protection against infective agents, chemical splashes and other particles. Manufactured of a soft and breathable material this coverall offers comfort and high-level protection. The coverall is equipped with hood, self-adhesive zipper, secure and convenient closures with bands on hoods, arms and legs. These coveralls can be applicable in various areas such as handling chemicals in dust and powder form, varnishing and protection against paint splashes, fiberglass production and processing, industrial cleaning and maintenance, automotive industry, grinding and polishing, mining and quarrying, demolition and refurbishment, working with asbestos, wood and metal processing, construction industry and pharmaceutical industry.

COLORS



SERVICES



SIZE	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	55	60	65	70	75	80	85	+/- 3 cm
B	188	190	192	194	196	198	200	+/- 3 cm
C	74	76	78	80	82	84	86	+/- 3 cm

SV- 113

DETAILS

Model: SV- 113 series
Fabric: Spunbond hotmelt polyethylen, PP + PE
Weight: 40-80 GSM
Cuff: Elastic
Color: White
Stitch: Overlock+Thermoband stripe
Other: Waist and neck straps for a secure fit
Optional: Static/Antistatic
Sterile/Non Sterile
Ultrasonic seam

For more informatin
please refer to the data sheets.




HI-VIS REFLECTIVE WORKWEAR COVERALL

Stay safe and visible in our High-Visibility Disposable Protective Coverall with Reflective Tape. Designed for both protection and visibility, this coverall is perfect for a wide range of applications, including construction, road work, industrial, and emergency response environments. Made from high-quality, breathable Microporous material, this coverall features reflective tape to enhance visibility in low-light situations, ensuring your safety while on the job.

COLORS



SERVICES




SIZE	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	55	60	65	70	75	80	85	+/- 3 cm
B	188	190	192	194	196	198	200	+/- 3 cm
C	74	76	78	80	82	84	86	+/- 3 cm

SV- 116

DETAILS

- Model: SV- 116 series
Fabric: SMS or Microporous PE Laminated
Weight: 40-90 GSM
Cuff: Elastic / Knitted cuff
Stitch: Overlock/ Ultrasonic stripe
Ergonomically shaped hood
Reflective tape
Anti-static
Optional: Sterile/Non Sterile

For more information,
please refer to the data sheets.



DISPOSABLE COVERALL (SMS)

These coveralls offer comfort and are durable, breathable and water-resistant. The fabric also silicone- and allergy-free. The coverall is equipped with hood, self-adhesive zipper, secure and convenient closures with bands on hoods, arms and legs. The welded stripes offer greater protection for liquids. These coveralls can be applicable in medical and industrial settings and conditions in which one works with dangerous biological and chemical substances.

SV- 188

DETAILS

Model: SV-188 series
Fabric: Spunbond+meltblown+spunbond
Weight: 20-80 GSM
Cuff: Elastic
Color: White
Stitch: Overlock

For more informatin
please refer to the data sheets.

DISPOSABLE COVERALL (SS)

These coveralls offer comfort and are durable, breathable and water-resistant. The fabric is also silicone- and allergy-free. The coverall is equipped with hood, self-adhesive zipper, secure and convenient closures with bands on hoods, arms and legs. These coveralls can be applicable in medical and industrial settings.


SV- 187

DETAILS

Model: SV-187 series
Fabric: Spunbond+Spunbond
Weight: 20-80 GSM
Cuff: Elastic
Color: White
Stitch: Overlock

For more informatin
please refer to the data sheets.

COLORS



SERVICES



Hoods & Bonnet

FOR SAFETY OF ALL IN HEALTHCARE SETTINGS

Disposable hoods ensure multiple protections for open wounds during surgeries. Additionally, disposable hoods avoid food to be contaminated by hair or dandruff in food industry or gastronomy.

 **PODiMA**

 **SUVICOM®**

 **rowimed®**

 **Raddax®**



SV- 49 | Round Elastic Bonnet Cap

Made of spunbond + polypropylene, soft and breathable fabric these caps be used by surgeons, nurses, physicians and other workers involved in patient care in hospitals. Equipped with elastic band around the head, these caps provide secure fit and offer complete hair coverage.

DETAILS

Fabric : SS , SMS, PP+PE
Weight : 20 - 60 gr

Double stitch
Anti_allergic.



SV- 50 | Operation Bonnet

Surgical bonnets are used to prevent contamination of the surgical, medical or hygiene sites by microorganisms originating from the hair of scalp of staff. These protective bonnets are applicable in medical, health, beauty and pharmaceutical settings.

DETAILS

Fabric : SS , SMS
Weight : 20 - 60 gr

Double stitch
Anti_allergic.



SV- 51 | Astro Hood

ASTRO Hoods are disposable surgical caps in astro-form, made of light material which protects from contamination. ASTRO hoods are applicable in high-contaminated medical settings, such as operating rooms. The size of an ASTRO Hood is individually adjustable.

DETAILS

Fabric : SS , SMS
Weight : 20 - 60 gr

Double stitch
Anti_allergic.



SV- 52 | Barret Hood

Breathable, dustproof disposable fleece caps which are made of environment- and skin-friendly material. Barret hoods offer effective protection against dust, fine dust, low-risk liquids and spraying. Ideal for hygiene, laboratory, industrial, medical settings.

DETAILS

Fabric : SS
Weight : 10 - 20 gr

Elastic
Anti_allergic.



SV- 53 | Surgical Hoods With Ties

Surgical hoods with ties are made of light-weight and comfortable fabric and provides, our hood provides complete head, neck and shoulder coverage. Adjustable back ties allow for a personalized fit.

DETAILS

Fabric : SS , SMS, PP+PE
Weight : 20 - 60 gr

Elastic
Anti_allergic.



SV- 54 | Headgear With Cape

Surgical hoods with ties are made of light-weight and comfortable fabric and provides, our hood provides complete head, neck and shoulder coverage. Adjustable back ties allow for a personalized fit.

DETAILS

Fabric : SS , SMS, PP+PE
Weight : 20 - 60 gr

Elastic
Anti_allergic.



SV- 55 | PE Cap

Disposable PE Shower Cap is made of superior quality, durable plastic, guaranteed not to tear during your bathing time. Enjoy your shower without fear of water getting through the transparent cap!

DETAILS

Fabric : PE
Weight : 9 - 20 Micron

Elastic
Anti_allergic.



Shoe Covers

FOR SAFETY OF ALL IN HEALTHCARE SETTINGS


Applicable in medical sector and food industry, shoe covers avoid contamination and offer protection both for ground (e.g. sandy shoes on damp ground) and shoes (e.g. protection from dust or paints during renovation works).

 **PODiMA**

 **SUVICOM®**

 **rowimed®**

 **Razdax®**



SV- 150 | Boot Cover

Long boot covers provide protection against infectious viruses. Long boot covers are tear-resistant non-slip premium quality with ergonomic construction that provides ease of use and better fit. At the same time, they offer a high level of breathability due to air permeability. They are made of materials repellent to water-based liquids and spray mist and are antistatic. Elasticity at the ankle allows one to put on the boot covers easily.

DETAILS

Fabric: PP+PE
Weight: 30-70 GSM
Stitch: Overlock
Optional: Ultrasonic seams, Non-slip velcro tape



SV- 151 | Thermo Band Striped Long Boot Cover

Long boot covers provide protection against infectious viruses. The heat-sealed tape on the seams prevents leakage and makes it the perfect accessory for our jumpsuit. Long boot covers are tear-resistant non-slip premium quality with ergonomic construction that provides ease of use and better fit. Elastic band and laces on the upper part provide extra coverage.

DETAILS

Fabric: PP+PE + Thermoband
Weight: 40-80 GSM
Stitch: Overlock
Optional: Ultrasonic seams, Non-slip velcro tape



SV- 152 | Short Overshoe

Overshoes complete our protective coveralls and provide additional protection against infectious viruses. Overshoes Laminated is microporous and its seams are hidden and do not permeate. They are made of soft durable, breathable and water-resistant fabric. They have non-slip soles, lace up and extra. They come with elastic band at the top for protection.

DETAILS

Fabric: PP+PE
Weight: 30-80 GSM
Stitch: Overlock
Optional: Ultrasonic seams, Non-slip velcro tape





SV- 60 | Shoe Cover Simple

Liquid-resistant overshoes are ideal for footwear hygiene and prevent potentially hazardous material from coming in contact with the bottom of a persons shoes. Thanks to the incorporated latex rubber, overshoes sit perfectly on feet and close tightly around ankles. With thicker materials, these overshoes are particularly easy to wear and do not tear quickly.

DETAILS

Fabric: PE Thickness: 9-40 Micron
 Stitch: Rope stitch technology is used to avoid shoe cover ladders
 Other: Reinforced double elastic band, liquid-proof, non-slip, tear-proof - Produced in automatic machines without being touched Compressed vacuum packaging



SV- 63 | SMS Long Boot Cover

This latex-free boot cover has breathable, fluid-resistant SMS material above the ankle, and a fluid-proof polypropylene coat below the ankle, providing maximum fluid protection. The heat-sealed seams reduce the chance of fluid penetration. The tread on the bottom helps provide excellent traction and increases stability. Knee High.

DETAILS

Fabric: SS or SMS
 Weight: 20-80 GSM
 Stitch: Overlock
 Optional: Ultrasonic seams, Non-slip velcro tape

Mask

FOR SAFETY OF ALL IN HEALTHCARE SETTINGS

Masks provide protection against harmful particles
thanks to their high filtering performance.

 **PODiMA**

 **SUVICOM®**

 **rowimed®**

 **Raddax®**



**HIGH FILTRATION
EFFICIENCY**



**HEALTHY
BREATHING**



**5-LAYER OF
PROTECTION**

FFP2 MASKS

FFP2 MASKS NR D

FFP2 Mask provides protection against harmful particles thanks to its high filtering performance. FFP2 mask with its water-repellent outer layer, is an effective personal protective equipment in protection from hazardous aerosols, gases and evaporations. It ensures convenient vision for eyeglass wearers and ease of use with its foldable shape.

COLORS



SERVICES



CE 2841 | **MADE IN
TÜRKİYE**

SV- 130 **SV- 131**
SV- 130-D **SV- 131-D**

DETAILS

Model: SV- 130-SV-131
SV-130-D - SV-131-D
Fabric: Spunbond + Meltblown
Sizes: Small / Medium / Standard
Weight: 150-180 GSM
Color: Black/White
Filtration: Efficiency \geq 95%

For more informatin
please refer to the data sheets.



SURGICAL MASKS


Surgical masks provide protection against infections and an array of pathogens and airborne particles. Without obstructing breath, the mask is made of light, soft and skin-friendly material and secure comfortable fit with integrated adjustable strip clip for nose bridge.

KIDS MASKS

Surgical masks provide protection against infections and an array of pathogens and airborne particles. Without obstructing breath, the mask is made of light, soft and skin-friendly material and secure comfortable fit with integrated adjustable strip clip for nose bridge.

2000R	2000R - K
<p>DETAILS SURGICAL MASKS</p> <p>Fabric: Spunbond + Meltblown Weight: 70-80 GSM Stitch: Ultrasonic seam Filtration: BFE > 95% Other: EU Type IIR I ASTM LVL 1</p> <p>Adjustable nose wire , Anti Allergic Air Permeable, Soft ear loops Splash resistant, Hygienic and hypoallergenic</p>	<p>DETAILS KIDS MASK</p> <p>Fabric: Spunbond + Meltblown Weight: 70-80 GSM Stitch: Ultrasonic seam Filtration: BFE > 95% Other: EU Type IIR I ASTM LVL 1</p> <p>Adjustable nose wire , Anti Allergic Air Permeable, Soft ear loops Splash resistant, Hygienic and hypoallergenic Optional Coloroud Furry design</p>
<p>For more informatin please refer to the data sheets.</p> <p>SERVICES</p> <div> </div>	

Disposable Bed Sheet



BED SHEET

Bed covers are ideal for any type of medical table, stretcher or patient bed. This includes hospital bed mattress cover for ambulance stretchers through to hospitality linens. They can be used in patient rooms in hospitals, clinics, or patient care facilities.

DETAILS

Fabric: SMS, SS
Weight: 35-60 GSM
Color: White/Green/Blue

COLORS



BED COVER

Made of non-woven fabric, disposable bed cover can reduce the risk of cross-infection effectively. Bed cover is soft, economic, lightweight, and breathable. Bed covers are also equipped with elastic bands on the two sides which is more convenient to use, and offer more protection. It is one of the best choices for patients' health in medical facilities.

DETAILS

Fabric: SMS, SS
Weight: 35-60 GSM
Color: White/Green/Blue

SERVICES



Disposable Towel & Sheet

HYGIENIC | ECO | ANTI BACTERIAL
ZERO BACTERIAL RISK

Absorbs water and dries body after showering.
Resists tearing when wet for more drying.
Great for travel, camping, long rides or daily use.
3 times absorbing more than cotton towels.

 **PODiMA**

 **SUVICOM®**

 **rowimed®**


 **Razdax®**

Disposable Towel



STANDARD
100

Biodegradable



SV- 630 | Hairdresser Towel **FOR WOMEN**

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for travel, camping, long rides, or daily use.

SIZE
40 x 80 cm
(16" x 31")

THICKNESS

STANDART ☒ THICK ☐
50gsm / m²

1 **DO NOT REUSE**
PER PACK

SV- 631 | Hairdresser Towel **FOR MEN**

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for travel, camping, long rides, or daily use.

SIZE
30x60 cm
(12" x 24")

THICKNESS

STANDART ☒ THICK ☐
50gsm / m²

1 **DO NOT REUSE**
PER PACK

SV- 632 | Feet Towel

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for travel, camping, long rides, or daily use.

SIZE
40 x 40 cm
(16" x 16")

THICKNESS

STANDART ☐ THICK ☒
70gsm / m²

1 **DO NOT REUSE**
PER PACK

SV- 633 | Hand Towel

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for travel, camping, long rides, or daily use.

SIZE
40 x 40 cm
(16" x 16")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 PER PACK **1** DO NOT REUSE

SV- 634 | Sports/Gym Towel

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for sport, travel, camping, long rides, or daily use.

SIZE
40 x 80 cm
(16" x 31")

THICKNESS
STANDART ☐ THICK ☒
70gsm / m²

1 PER PACK **1** DO NOT REUSE

SV- 635 | Bath Towel

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for travel, camping, long rides, or daily use.

SIZE
60 x 120 cm
(24" x 47")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 PER PACK **1** DO NOT REUSE

SV- 636 | Bath Towel

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for Sport , travel, camping, long rides, or daily use.

SIZE
80x140 cm
(31" x 55")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 PER PACK **1** DO NOT REUSE

SV- 637 | Pet Towel FOR FRIENDS

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for Pet, travel, camping, long rides, or daily use.



SIZE
80 x 140 cm
(31" x 55")

THICKNESS
STANDART ☐ THICK ☒
70gsm / m²

1 PER PACK **1** DO NOT REUSE

Disposable Towel Set

SV- 638 | Towel Set

Set of 3 (Hair, bath, foot)

Absorbs water and dries body after showering Resists tearing when wet for more drying Great for travel, camping, long rides, or daily use.

SIZE / HAIR
40 x 80 cm
(16" x 31")

SIZE / BATH
80 x 140 cm
(31" x 55")

SIZE / FOOT
40 x 45 cm
(16" x 18")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 PER PACK **1** DO NOT REUSE



SV- 639 | Towel Set

Set of 2 (Hair, bath)

Absorbs water and dries body after showering Resists tearing when wet for more drying Great for travel, camping, long rides, or daily use.

SIZE / HAIR
40 x 80 cm
(16" x 31")

SIZE / BATH
80 x 140 cm
(31" x 55")

THICKNESS
STANDART ☐ THICK ☒
70gsm / m²

1 PER PACK **1** DO NOT REUSE

Disposable Bed Set

Biodegradable



STANDARD
100

Personal use disposable bed set is hygienic and ultra light. It doesn't contain any synthetic material and doesn't cause sweat. It holds zero bacteria risk.

SV- 640 Single Bed Set 1 Bed sheet, 1 pillowcase, 1 duvet cover		SV- 641 Double Bed Set 1 Bed sheet, 2 pillowcase, 1 duvet cover	
Disposable duvet cover set is hygienic and ultra-light. Travel, Hospital, hotel, camping, caravan, yacht, Personal use...		Disposable duvet cover set is hygienic and ultra-light. Travel, Hospital, hotel, camping, caravan, yacht, Personal use...	
BED SHEET 160x 200 cm (64" x 80")		BED SHEET 200x 240 cm (80" x 96")	
PILLOW CASE 50x 70 cm (20" x 28")		PILLOW CASE 50x 70 cm (20" x 28")	
THICKNESS STANDART <input checked="" type="checkbox"/> THICK <input type="checkbox"/> 50gsm / m ²		THICKNESS STANDART <input checked="" type="checkbox"/> THICK <input type="checkbox"/> 50gsm / m ²	
1 SET PACK		1 SET PACK	
DO NOT REUSE		DO NOT REUSE	
DUVET COVER 160x 200 cm (64" x 80")		DUVET COVER 160x 200 cm (80" x 96")	

Disposable Travel Set

Personal use disposable bed set is hygienic and ultra light. It doesn't contain any synthetic material and doesn't cause sweat. It holds zero bacteria risk.

SV- 642 Single Travel Set 1 Bath towel, 1 Hair towel, 1 Foot towel 1 Bed Sheet, 1 Pillow Case		BATH TOWEL 80x 140 cm (31" x 55")		FOOT TOWEL 40x 45 cm (16" x 18")	
Disposable duvet cover set is hygienic and ultra-light. Travel, Hospital, hotel, camping, caravan, yacht, Personal use...		HAIR TOWEL 40x 80 cm (16" x 31")		PILLOW CASE 50x 70 cm (20" x 28")	
THICKNESS STANDART THICK <input checked="" type="checkbox"/> <input type="checkbox"/> 50gsm / m ²		DUVET COVER 160x 200 cm (63" x 79")			
1 SET PACK					
DO NOT REUSE					

SV- 643 Double Travel Set 2 Bath towel, 2 Hair towel, 2 Foot towel 1 Bed Sheet, 2 Pillow Case		BATH TOWEL 80x 140 cm (31" x 55")		FOOT TOWEL 40x 45 cm (16" x 18")	
Disposable duvet cover set is hygienic and ultra-light. Travel, Hospital, hotel, camping, caravan, yacht, Personal use...		HAIR TOWEL 40x 80 cm (16" x 31")		PILLOW CASE 50x 70 cm (20" x 28")	
THICKNESS STANDART THICK <input checked="" type="checkbox"/> <input type="checkbox"/> 50gsm / m ²		DUVET COVER 160x 200 cm (63" x 79")			
1 SET PACK					
DO NOT REUSE					



SV- 644 | Double Duvet Case

Disposable duvet cover set is hygienic and ultra-light. Travel, hospital, hotel, camping, caravan, yacht, Personal use...

SIZE
200 x 240 cm
(79" x 96")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 
PER PACK DO NOT REUSE

SV- 645 | Single Duvet Case

Disposable duvet cover set is hygienic and ultra-light. Travel, hospital, hotel, camping, caravan, yacht, Personal use...

SIZE
160 x 200 cm
(63" x 79")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 
PER PACK DO NOT REUSE

SV- 646 | Double Bed Sheet

Disposable duvet cover set is hygienic and ultra-light. Travel, hospital, hotel, camping, caravan, yacht, Personal use...

SIZE
200 x 240 cm
(79" x 96")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 
PER PACK DO NOT REUSE

SV- 647 | Single Bed Sheet

Disposable duvet cover set is hygienic and ultra-light. Travel, hospital, hotel, camping, caravan, yacht, Personal use...

SIZE
160 x 200 cm
(63" x 79")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 
PER PACK DO NOT REUSE

SV- 648 | Single Pillow Case

Disposable duvet cover set is hygienic and ultra-light. Travel, hospital, hotel, camping, caravan, yacht, Personal use...

SIZE
50 x 80 cm
(20" x 31")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 
PER PACK DO NOT REUSE

Disposable Massage Sheet & Sun longer Cover

Biodegradable

SV- 649 | Spa Massage Sheet

Disposable duvet cover set is hygienic and ultra-light. Hospital, hotel, yacht, Professional use...

SIZE
80 x 200 cm
(31" x 79")

THICKNESS

STANDART ☒ THICK ☐
50gsm / m²

1
PER PACK DO NOT REUSE

SV- 650 | Spa Massage Sheet (hole)

Disposable duvet cover set is hygienic and ultra-light. Hospital, hotel, yacht, Professional use...

SIZE
80 x 200 cm
(31" x 79")

THICKNESS

STANDART ☒ THICK ☐
50gsm / m²

1
PER PACK DO NOT REUSE

SV- 651 | Sun Longer Cover

Disposable duvet cover set is hygienic and ultra-light. Hospital, hotel, yacht, Professional use...

SIZE
80 x 220 cm
(31" x 87")

THICKNESS

STANDART ☒ THICK ☐
50gsm / m²

1
PER PACK DO NOT REUSE

SV- 652 | Sun Longer Cover with pouch

Disposable duvet cover set is hygienic and ultra-light. Beach, hotel, yacht, Professional use...

SIZE
80 x 200 cm
(31" x 79")

THICKNESS

STANDART ☒ THICK ☐
50gsm / m²

1
PER PACK DO NOT REUSE

Our Certificates

Certificates Modul B



Certificates Modul C2



Certificates




Test Reports



Production Process



■ Raw Material



■ Tailoring



■ Sewing



■ Pressure Plastic Molding



■ Packaging



■ Delivery

How is production carried out

Production process and management are conducted based on the internationally accepted Strict hygiene standards. The products are duly packaged and delivered to the end users on time following the guiding hygiene principles.



Fast delivery

Fast delivery from our Europe central warehouse with high capacity Turkey based production.

EU certification

In order to be able to offer the best protection and the greatest security in everyday working life, we only supply certified products.

Top products

We only carry tested, certified and high quality products.

Fair prices

Thanks to strong partnerships and optimized business processes, we can deliver on our products at fair prices.

Declaration On Product Standards & Services



CATEGORY 3

Liquid-tight protective clothing: Full body or full protection suits against liquid chemicals



EN 13688

General requirement for protective clothing.
EN ISO 13688:2013



EN 14126

Protective clothing against infectious agents.
EN 14126: 2003 + AC:2004



EN 14605 TYPE PB 3-B

Liquid-tight protective clothing against liquid chemicals.
EN 14605:2005 + A1:2009 TYP PB 3-B



EN 13034 TYPE PB 6-B

Protective clothing against liquid chemicals with limited protection.
EN 13034:2005 + A1:2009 TYP PB 6-B



EN 14605 TYPE PB 4-B

Spray-tight protective clothing against liquid chemicals.
EN 14605:2005+A1:2009



EN 13982 TYPE 5-B

Protective clothing against airborne solid particles
EN ISO 13982-2:2004 TYP 5



EN 1149

Protective clothing with electrostatic properties
EN 1149-1:2006

EN 13795

Surgical drapes, gowns and clean air suits, used as medical devices for patients, clinical staff and equipment.



EN 1073-2:2002

Protective clothing against radioactive contamination-part2

FFP2 NR D

Passes EN 149:2001 optional Dolomite (D) clogging test for more comfortable breathing.

Services



HOSPITAL



INDUSTRY



CORONA PROTECTION



EMERGENCY



HYGIENE



LABORATORY



About Podima

Our products are produced in accordance with our production infrastructure, expert staff and strict hygiene standards in Türkiye. Accordingly, our products are used and demanded safely in Europe and all over the world. Our quality, production and distribution performance has been appreciated by all our customers, especially during the pandemic period.

Our **SUVICOM** brand is primarily used in Germany and our other brand, **ROWIMED**, is used in Italy with its high quality and with confidence.

Our Philosophy Our goal is to supply quality medical PPEs and supplies reliably and continuously at fair prices. We are committed to ensure a gap-free supply chain system, and achieve highest customer satisfaction and sustainable partnerships.

Aiming highest supply standards, quality is key for us, which is ensured through constant controls and quality checks. We regularly do site visits to our manufacturers in order to personally convince ourselves of the production facilities, employees, materials and business processes. Only products duly certified are included in our portfolio.

Imprint / Contact



Head Office
Podima Medikal ve
Tekstil San. Tic. Ltd. Şti.
Şerifali Mh, Mevdudi Sk.
No:53/1, 34775 Ümraniye
İstanbul / Türkiye

Phone: +90 216 784 07 93
+90 216 314 93 16
Email: info@podimamedikal.com

■ Factory:
Siirt Yolu Cad. 72100
Merkez Batman / Türkiye

■ Branch Office
Suivicom Media GmbH
Lilienthalstr 14
41515, Grevenbroich,
Germany

Phone: +49 221 292033 31
Email: info@suivicom.de
Web: www.suivicom-medical.de