



EU-MDR

Protective Clothing

Your reliable partner for personal protective equipment, including coveralls, protective gowns, hoods, shoe covers, nitrile gloves and protective masks.

CE 0425 **OEKO TEX®** STANDARD 100 **CE** 2841



Gowns

Disposable visitor, examination and surgical gowns for medical sector to protect health professionals, patients and visitors against viruses and germs. Also applicable in food industry, manufacturing, laboratories and pharmacies.



Coveralls

Disposable coveralls ensure protection against infective agents, chemical splashes and other particles. Manufactured of a soft and breathable material coveralls offer comfort and full-body protection.



Towel & Sheets



Hoods

Disposable hoods ensure multiple protections for open wounds during surgeries. Additionally, disposable hoods avoid food to be contaminated by hair or dandruff in food industry or gastronomy.




Shoe covers

Applicable in medical sector and food industry, shoe covers avoid contamination and offer protection both for ground (e.g. sandy shoes on damp ground) and shoes (e.g. protection from dust or paints during renovation works).

CONTENTS


NEUROSURGERY DRAPES.....	4-5
GENERAL SURGERY PACK.....	6-8
GYNECOLOGY PACK.....	9-11
E.N.T. PACK.....	12-13
CARDIOLOGY PACK.....	14-15
CARDIOVASCULAR PACK.....	16
OPHTALMIC PACK.....	17
ORTHOPEDICS PACK.....	18-20
UROLOGY PACK.....	21-22
OTHER PACK & DRAPS.....	23-29
DISPOSABLE GOWNS.....	30-35
DISPOSABLE COVERALLS.....	36-37
DISPOSABLE FFP2 - FFP3 MASK & SURGICAL MASK.....	38
DISPOSABLE BONNET & SHOE COVER.....	39
DISPOSABLE TOWELS & SHEET.....	40-44
OUR CERTIFICATES.....	45
PRODUCTION PROCESS.....	46
STANDARDS.....	47

CERTIFICATIONS AND STANDARDS



ISO 13485:2016 Certified	ISO 9001-2015 Certified	ISO 45001-2018 Certified
ISO 10002-2018 Certified	ISO 14001-2015 Certified	


NEUROSURGERY / BEYİN CERRAHİ



REF SV-1010

NO	QTY	NAME	SIZE
1	3	Towel / Havlu	40x40
2	1	Reinforced Surgical Gown / Korumalı cerrahi önlük	M
3	2	Reinforced Surgical Gown / Korumalı cerrahi önlük	XL
4	4	Side Banded Drape / Yan Bantlı Örtü	75x90
5	1	Craniotomy Drape / Kraniotomi Örtüsü	240x320
6	1	Basic Drape / Düz Örtü	150x250
7	10	Gauze Swab / Gaz Spanç	7,5x7,5
8	1	Skin Marker And Ruler / Cetvel ve Cilt Kalemi	ST
9	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
10	2	Back Table Cover / Alet Masa Örtüsü	150x200
11	1	Sterilization Wrap / Sterilizasyon Boşçası	100x100


CARPAL TUNNEL PACK / KARPAL TUNEL SETİ



REF SV-1011

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı cerrahi önlük	L
3	1	Anesthesia Drape / Anestezi Örtüsü	150x240
4	1	Basic Drape / Düz Örtü	160x180
5	1	Carpal Tunnel Drape / Karpal Tunel Örtüsü	150x150
6	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
7	2	Back Table Cover / Alet Masa Örtüsü	150x200
8	1	Sterilization Wrap / Sterilizasyon Boşçası	100x100


HEAD PACK / KAFA SETİ



REF SV-1012

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	1	Head (Turban) Drape / Baş Örtüsü (Türban)	100x120
3	1	Foot Drape / Bantlı Ayak Örtüsü	150x180
4	1	Fenestrated Drape (10cm with adhesive) / Delikli Bantlı Örtü / (10 cm)	100x100
5	1	Back Table Cover / Alet Masa Örtüsü	150x200
6	1	Sterilization Wrap / Sterilizasyon Boşçası	50x50


LAMINECTOMY PACK / LAMİNOKTEMİ SETİ



REF SV-1013

NO	QTY	NAME	SIZE
1	3	Towel / Havlu	40x40
2	1	Standard Surgical Gown / Standart Cerrahi Önlük	M
3	2	Standard Surgical Gown / Standart Cerrahi Önlük	XL
4	4	Side Banded Drape / Yan Bantlı Örtü	75x90
5	1	Spinal Vertebra Drape With Incise Film / Spinal Vertebra Örtüsü (İnsizyon Film)	150x300
6	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
7	1	Back Table Cover / Alet Masa Örtüsü	150x200
8	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


CRANIOTOMY PACK / KRANIOTOMİ SETİ



REF SV-1014

NO	QTY	NAME	SIZE
1	1	Craniotomy Drape / Kraniyotomi Örtüsü	240x320
2	1	Back Table Cover / Alet Masa Örtüsü	150x200
3	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
4	1	Basic Drape / Düz Örtü	150x200
5	4	Side Banded Drape / Yan Bantlı Örtü	50x50
6	2	Op-Tape / Yapıştırma Bandı	10x50
7	1	Crepe Paper / Krep Kağıdı	100x100
8	2	Towel / Havlu	40x40


SHUNT PACK / SHUNT SETİ



REF SV-1015

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	XL
3	4	Side Banded Drape / Yan Bantlı Örtü	75x90
4	1	Anesthesia Drape / Anestezi Örtüsü	150x240
5	1	Shunt Drape / Shunt Örtüsü	150x300
6	1	Incise Film / Drep (İnsizyonlu Film)	25x75
7	1	Basic Drape / Düz Örtü	150x250
8	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
9	1	Back Table Cover / Alet Masa Örtüsü	150x200
10	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


GENERAL SURGERY PACK / GENEL CERRAHİ SETİ



REF SV-1016

NO	QTY	NAME	SIZE
1	4	Towel / Havlu	40x40
2	3	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	1	Anesthesia Drape / Anestezi Örtüsü	150x240
4	2	Side Banded Drape / Yan Bantlı Örtü	75x90
5	1	Foot Drape / Ayak Örtüsü	150x180
6	1	Op Tape / Yapıştırma Bandı	10x50
7	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
8	2	Back Table Cover / Alet Masa Örtüsü	150x200
9	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


GENERAL SURGERY PACK (STANDARD) / GENEL CERRAHİ SETİ (STANDART)



REF SV-1017

NO	QTY	NAME	SIZE
1	4	Towel / Havlu	40x40
2	1	Anesthesia Drape / Anestezi Örtüsü	150x240
3	2	Side Banded Drape / Yan Bantlı Örtü	75x90
4	1	Foot Drape / Ayak Örtüsü Bantlı	150x180
5	1	Op-Tape / Yapıştırma Bandı	10x50
6	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
7	1	Back Table Cover / Alet Masa Örtüsü	150x200
8	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


GENERAL SURGERY PACK (BASIC) / GENEL CERRAHİ SETİ (TEMEL)



REF SV-1018


NO	QTY	NAME	SIZE
1	4	Towel / Havlu	40x40
2	1	Anesthesia Drape / Anestezi Örtüsü	150x240
3	2	Side Banded Drape / Yan Bantlı Örtü	75x90
4	1	Foot Drape / Ayak Örtüsü Bantlı	150x180
5	1	Op-Tape / Yapıştırma Bandı	10x50
6	1	Back Table Cover / Alet Masa Örtüsü	100x150
7	1	Sterilization Wrap / Sterilizasyon Bohçası	75x75

OPERATION PACK / CERRAHİ MÜDAHALE SETİ



REF SV-1019


NO	QTY	NAME	SIZE
1	1	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
2	1	Fenestrated Drape / Delikli Örtü	75x75
3	1	Back Table Cover / Alet Masa Örtüsü	100x120

LAPAROSCOPIC PELVISCOPY PACK /
LAPAROSKOPİ PELVİSKOPİ SETİ

REF SV-1020

NO	QTY	NAME	SIZE
1	4	Towel / Havlu	40x40
2	1	Tube Holder [Velcro] / Tüp Tutucu [Bantlı Cırt]	2.5x30
3	1	Mobile Pouch PE / Seyyar Poşet PE	28x33
4	1	Laparoscopic Pelviscopy Abdominal Drape Litotomic / Laparoskopî Pelviskopî Örtü Litotomik	180/240x260
5	1	Back Table Cover / Alet Masa Örtüsü	150x190
6	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
7	1	Op-Tape / Yapıştırma Bandı	10x50


LAPARATOMY PACK / LAPAROTOMİ SETİ



REF SV-1021

NO	QTY	NAME	SIZE
1	6	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	20	Gauze Swab / Gaz Spanç	7.5X7.5
4	1	Laparotomy Drape With Incise Film(25x30) / Laparatomi Örtüsü (25x30) İnsizyonlu	240X360
5	1	Mayo Stand Cover / Mayo Masa Örtüsü	78X145
6	1	Back Table Cover / Alet Masa Örtüsü	150X200
7	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


LAPAROTOMY PACK WITH POUCH / LAPAROTOMİ SETİ POŞLU



REF SV-1022

NO	QTY	NAME	SIZE
1	6	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	20	Gauze Swab / Gaz Spanç	7.5X7.5
4	1	Laparotomy Drape With Incise Film And Pouch [25x30] / Laparotomi Örtüsü (25x30) İnsizyonlu Poşlu	240X360
5	1	Mayo Stand Cover / Mayo Masa Örtüsü	78X145
6	1	Back Table Cover / Alet Masa Örtüsü	150X200
7	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


LAPAROSCOPIC ABDOMINAL PERINEAL PACK / LAPAROSKOPİ ABDOMİNAL PERİNEAL SETİ



REF SV-1023

NO	QTY	NAME	SIZE
1	4	Towel / Havlu	40x40
2	1	Side Banded Drape / Yan Bantlı Örtü	75x75
3	1	Laparoscopic Abdominal Perineal Drape [25x30 and 12x16 Fenestration] Litotomic/ Laparoskopik Abdominal Perineal Örtü(25x30 ve 12x16 Delik Bant)	180/240x260
4	1	Back Table Cover / Alet Masa Örtüsü	150x190
5	1	Op-Tape / Yapıştırma Bandı	10x50
6	1	Tube Holder (Velcro) / Tüp Tutucu (Bantlı Cırt)	2.5x30


EMERGENCY DELIVERY PACK / ACİL DOĞUM SETİ



REF SV-1024

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	1	Standard Surgical Gown / Standart Cerrahi Önlük	L
3	1	Under Buttocks Drape / Hasta Alt Örtüsü	75x150
4	1	Mask / Maske	ST
5	1	Baby Drape / Bebek Örtüsü	100x100
6	1	Abdominal Drape / Abdominal Örtü	100x120
7	2	Gloves (Pair) / Eldiven (Çift)	7 Nr.
8	10	Gauze Swab / Gaz Spanç	7.5x7.5
9	1	Umbilical Clamp / Göbek Klempi	ST
10	2	Maternity Pad / Doğum Pedi	L
11	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100

E.T PACK / E.T SETİ



REF SV-1025

NO	QTY	NAME	SIZE
1	1	OPU Drape / O.P.U Örtüsü (Delikli)	130x150
2	1	Basic Drape / Düz Örtü	80x90
3	10	Gauze Swab / Gaz Spanç	7.5x7.5
4	1	Sterilization Wrap / Sterilizasyon Bohçası	75x75


GYNECOLOGY O.P.U. PACK / JİNEKOLOJİ O.P.U. SETİ



REF SV-1026

NO	QTY	NAME	SIZE
1	1	Towel / Havlu	40x40
2	2	Standard Surgical Gown / Standart Cerrahi Önlük	L
3	2	Lithotomic Leg Cover / Litotomik Bacak kılıfı	55x120
4	1	OPU Drape With Fenestration / O.P.U Örtüsü (Delikli)	135x150
5	1	Plain Drape / Düz Örtü	80x90
6	1	Transparent PE Adhesive Drape / Yapışkan Bantlı Şeffaf Örtü PE	60x60
7	10	Gauze Swab / Gaz Spanç	7.5x7.5
8	1	Back Table Cover / Alet Masa Örtüsü	100x130
9	1	Sterilization Wrap / Sterilizasyon Bohçası	75x75


GYNECOLOGY PACK / JİNEKOLOJİ SETİ



REF SV-1027

NO	QTY	NAME	SIZE
1	6	Towel / Havlu	40x40
2	1	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	2	Leg Cover / Litotomik Bacak Kılıfı	55x120
4	1	Back Table Cover / Alet Masa Örtüsü	150x200
5	1	Abdominal Drape / Abdominal Örtü	100x120
6	1	Gynaecology Split Drape / Jinekoloji Örtüsü	75x120
7	1	Basic Drape / Düz Örtü	80x90
8	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


DELIVERY PACK / KADIN DOĞUM SETİ



REF SV-1028

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	1	Under Buttocks Drape / Hasta Alt Örtüsü	100x150
4	1	Abdominal Drape / Abdominal Örtü	75X100
5	1	Baby Drape / Bebek Örtüsü	100x100
6	2	Lithotomic Leg Cover / Litotomik Bacak kılıfı	55X120
7	1	Back Table Cover / Alet Masa Örtüsü	100X130
8	1	Sterilization Wrap / Sterilizasyon Bohçası	100X100


DELIVERY PACK (WITH POUCH) / POŞLU DOĞUM SETİ



REF SV-1029

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	1	Under Buttocks Drape With Pouch / Hasta Alt Örtüsü Poşlu	100x130
3	1	Abdominal Drape / Abdominal Örtü	75x100
4	1	Baby Drape / Bebek Örtüsü	100x100
5	2	Lithotomic Leg Cover / Litotomik Bacak kılıfı	55X120
6	1	Back Table Cover / Alet Masa Örtüsü	100x150
7	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


CESAREAN SECTION PACK / L.T.T. SEZARYEN SETİ



REF SV-1030

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	1	Cesarean Section Drape With Pouch / Sezaryen Örtüsü L.T.T	200x300
4	1	Baby Drape / Bebek Örtüsü	100x100
5	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
6	2	Umbilical Clamp / Göbek Klempİ	ST
7	1	Back Table Cover / Alet Masa Örtüsü	150x200
8	1	Sterilization Wrap / Sterilizasyon Bohçası	100X100


LITHOTOMY GYNECOLOGY PACK / LİTOTOMİ JİNEKOLOJİ SETİ



REF SV-1031

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Standard Surgical Gown / Standart Cerrahi Önlük	L
3	1	Lithotomic Gynaecology Drape / Litotomik Jinekoloji Örtüsü	180x240x240
4	1	Back Table Cover / Alet Masa Örtüsü	150x200
5	1	Tube Holder (Velcro) / Tüp Tutucu (Bantlı Cırt)	2.5x30
6	1	Sterilization Wrap / Sterilizasyon Bohçası	75x75


E.N.T. PACK (THROAT) / KBB SETİ (BOĞAZ)



REF SV-1032

NO	QTY	NAME	SIZE
1	1	Towel / Havlu	40x40
2	1	Reinforced Surgical Gown / Korumalı cerrahi önlük	L
3	1	E.N.T. Drape / KBB Örtüsü	150x240
4	1	Fenestrated Drape / Delikli Bantlı Örtü	100x100
5	1	Basic Drape / Düz Örtü	150x200
6	1	Back Table Cover / Alet Masa Örtüsü	100x150
7	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


E.N.T. PACK (NOSE) / KBB SETİ (BURUN)



REF SV-1033

NO	QTY	NAME	SIZE
1	1	Towel / Havlu	40x40
2	1	Reinforced Surgical Gown / Korumalı cerrahi önlük	L
3	1	E.N.T. Drape / KBB Örtüsü	150x240
4	1	Basic Drape / Düz Örtü	150x200
5	1	Back Table Cover / Alet Masa Örtüsü	100x150
6	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


TONSIL ADENOID PACK / TONSİL ADENOİT SETİ



REF SV-1034


NO	QTY	NAME	SIZE
1	1	Towel / Havlu	40x40
2	1	Standard Surgical Gown / Standart Cerrahi Önlük	L
3	1	Fenestrated Drape / Delikli Bantlı Örtü	100x100
4	2	Side Banded Drape / Yan Bantlı Örtü	150x240
5	1	Basic Drape / Düz Örtü	100x120
6	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
7	1	Back Table Cover / Alet Masa Örtüsü	150x200

THYROID PACK / TROİD SETİ




REF SV-1035

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı cerrahi önlük	L
3	1	Op Tape / Yapıştırma Bandı	10x30
4	1	Thyroid Drape / Troid Örtüsü	240x320
5	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
6	1	Back Table Cover / Alet Masa Örtüsü	150x200
7	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


**ANGIOGRAPY PACK-1/
ANJİYOGRAFİ SETİ-1**

REF SV-1036

NO	QTY	NAME	SIZE
1	2	Towel/Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı cerrahi önlük	L
3	2	Femoral Angiography Drape / Femoral Anjiyo Örtüsü	150x300
4	1	Basic Drape / Düz Örtü	100x200
5	1	Gallipot / Sıvı Kabı	250ml
6	1	Fluoroscopy Cover / Fluoroskopi Kılıfı	90Q
7	10	Gauze Swab / Gaz Spanç	7.5x7.5
8	1	Betadine Sponge / Boyama Sünger	ST
9	1	Back Table Cover / Alet Masa Örtüsü	150x150
10	1	Sterilization Wrap / Sterilizasyon Boşçası	100x100

**ANGIOGRAPY PACK-2 /
ANJİYOGRAFİ SETİ-2**

REF SV-1037

NO	QTY	NAME	SIZE
1	2	Towel/Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı cerrahi önlük	L
3	1	Plain Drape / Düz Örtü	100x200
4	1	Angiography Drape With Transparent Panel / Panelli Anjiyo Örtüsü	180x300
5	1	Fluoroscopy Cover / Fluoroskopi Kılıfı	90Q
6	1	Gallipot / Sıvı Kabı	250 ml
7	1	Gallipot / Sıvı Kabı	500 ml
8	10	Gauze Swab / Gaz Spanç	7.5x7.5
9	1	Betadine Sponge / Boyama Sünger	ST
10	2	Syringe / Enjektör	10 CC
11	1	Function Needle / Fonksiyon İğnesi	18
12	1	Opaque Custom Set (3 Taps) / Opak Set (3 Musluklu)	ST
13	1	6F Intra Ducer / 6F Intra Ducer	ST
14	1	Guidewire / Kılavuz Tel	ST
15	1	Surgical Blade / Bisturi Ucu	11
16	1	Tray / Tepsi	2500 ml
17	1	Back Table Cover / Alet Masa Örtüsü	150x150
18	1	Sterilization Wrap / Sterilizasyon Boşçası	100x100


ANGIOGRAPY PACK-3/ ANJİYOGRAFİ SETİ-3



REF SV-1038

NO	QTY	NAME	SIZE
1	5	Towel/Havlu	40x40
2	1	Femoral- Radial Angiography Drape With Transparent Panels (120x300 65 cm Panels, Femoral 12 cm With Incision 7 cm Radial 7x12 cm Oval Fenestration 120x100 cm Reinforcement) / "Femoral-Radyal Anjiyografi Örtüsü 4 Açıklık (Femoral 12 cm, Radyal 7x12 Oval İnsizyon) Çift Taraf Panelli	220x300
3	1	Basic Drape / Düz Örtü	100x200
4	1	Gallipot / Sıvı Kabı	500ml
5	1	Gallipot / Sıvı Kabı	250ml
6	1	Fluoroscopy Cover / Fluoroskopi Kılıfı	90Q
7	10	Gauze Swab / Gaz Spanç	7.5x7.5
8	1	Betadine Sponge / Boyama Sünger	ST
9	1	Back Table Cover /Alet Masa Örtüsü	150x150
10	1	Sterilization Wrap/ Sterilizasyon Bohçası	100x100

BIOPSY PACK / BİOPSİ SETİ



REF SV-1039

NO	QTY	NAME	SIZE
1	1	Fenestrated Drape / Delikli Örtü	80x100
2	1	Kidney Tray / Böbrek Küveti	ST
3	10	Gauze Swab / Gaz Spanç	7,5x7,5
4	1	Back Table Cover /Alet Masa Örtüsü	50x60
5	1	Sterilization Wrap/ Sterilizasyon Bohçası	50x50


BY-PASS PACK / BY-PASS SETİ



REF SV-1040

NO	QTY	NAME	SIZE
1	5	Towel / Havlu	40x40
2	1	Reinforced Surgical Gown / Korumalı Cerrahi önlük	M
3	2	Reinforced Surgical Gown / Korumalı Cerrahi önlük	L
4	1	Cardiovascular Drape / KVC Örtüsü	240x200x340
5	1	Side Banded Drape / Yan Bantlı Örtü	75x90
6	1	"U" Split Perineum Drape / "U" Yanklı Perine Örtüsü	150x260
7	2	Feet Socks / Patik	ST
8	2	Mobile Pouch PE / Seyyar Poşet PE	28x33
9	2	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
10	1	Back Table Cover / Alet Masa Örtüsü	150x200
11	2	Op-Tape / Yapışkan Bant	10x50
12	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


HEART VALVE SURGERY PACK / KALP KAPAK SETİ



REF SV-1041

NO	QTY	NAME	SIZE
1	5	Towel / Havlu	40x40
2	1	Reinforced Surgical Gown / Korumalı Cerrahi önlük	M
3	2	Reinforced Surgical Gown / Korumalı Cerrahi önlük	L
4	1	Cardiovascular Drape / KVC Örtüsü	240x200x340
5	1	Side Banded Drape / Yan Bantlı Örtü	75x90
6	2	Mobile Pouch PE / Seyyar Poşet PE	28x33
7	2	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
8	1	Back Table Cover / Alet Masa Örtüsü	150x200
9	2	Op-Tape / Yapışkan Bant	10x50
10	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


CONGENITAL CARDIOVASCULAR PACK / KONGENİTAL KARDİYOVASKÜLER SETİ



REF SV-1042

NO	QTY	NAME	SIZE
1	5	Towel / Havlu	40x40
2	1	Reinforced Surgical Gown / Korumalı Cerrahi önlük	M
3	2	Reinforced Surgical Gown / Korumalı Cerrahi önlük	L
4	1	Cardiovascular Drape / KVC Örtüsü	240x340
5	1	Side Banded Drape / Yan Bantlı Örtü	75x90
6	2	Mobile Pouch PE / Seyyar Poşet PE	28x33
7	2	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
8	1	Back Table Cover / Alet Masa Örtüsü	150x200
9	2	Op-Tape / Yapışkan Bant	10x50
10	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


OPHTALMIC PACK WITH SINGLE POUCH / OFTALMİK SET TEK POŞLU



REF SV-1043

NO	QTY	NAME	SIZE
1	2	Towel /Havlu	40x40
2	2	Standard Surgical Gown / Standart Cerrahi Önlük	L
3	1	Opthalmic Drape With Single Pouch / Tek Poşlu Oftalmik Örtü	100X130
4	1	Basic Drape / Düz Örtü	150x180
5	1	Back Table Cover / Alet Masa Örtüsü	100x130
6	1	Sterilization Wrap / Sterilizasyon Bohçası	75x75


OPHTALMIC PACK WITH DOUBLE POUCH / OFTALMİK SET ÇİFT POŞLU



REF SV-1044

NO	QTY	NAME	SIZE
1	2	Towel /Havlu	40x40
2	2	Standard Surgical Gown / Standart Cerrahi Önlük	L
3	1	Opthalmic Drape With Double Pouch / Çift Poşlu Oftalmik Örtü	150x200
4	1	Basic Drape / Düz Örtü	150x180
5	1	Back Table Cover / Alet Masa Örtüsü	100x130
6	1	Sterilization Wrap / Sterilizasyon Bohçası	75x75

ARTHROSCOPY PACK WITHOUT POUCH / ARTROSKOPİ SETİ POŞSUZ



REF SV-1045

NO	QTY	NAME	SIZE
1	2	Towel /Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	1	Sub Extremity Drape / Alt Extremitte Örtüsü	200x310
4	1	Leg Cover / Bacak Kılıfı	40x75
5	1	Elastic Bandage / Elastik Bandaj	10x150
6	1	Op Tape / Yapışkan Bant	10x50
7	1	Back Table Cover / Alet Masa Örtüsü	150x200
8	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


ARTHROSCOPY PACK WITH POUCH / ARTROSKOPİ SETİ POŞLU



REF SV-1046

NO	QTY	NAME	SIZE
1	2	Towel /Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	1	Arthroscopy Drape With Liquid Pouch / Artroskopî Örtüsü Likit Toplama Poşlu	200X300
4	1	Leg Cover / Bacak Kılıfı	40X75
5	1	Elastic Bandage / Elastik Bandaj	10X150
6	1	Skin Marker And Ruler / Cetvel ve Cilt Kalemî	ST
7	1	Tube Cover / Kamera Kılıfı	14X250
8	1	Op Tape / Yapışkan Bant	10x50
9	1	Back Table Cover / Alet Masa Örtüsü	150x200
10	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


HIP PACK / HİP SETİ



REF SV-1047

NO	QTY	NAME	SIZE
1	2	Towel /Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	2	Side Banded Drape / Yan Bantlı Örtü	100X120
4	1	"U" Split Drape PE / "U" Yanklı Örtü PE	150X200
5	1	Anesthesia Drape / Anestezi Örtüsü	150X240
6	1	"U" Split Drape / "U" Yanklı Örtü	150X300
7	1	Foot Drape / Ayak Örtüsü	160X180
8	2	X-Ray Cassette Bag / Rontgen Kaset Torbası	60X80
9	1	Back Table Cover / Alet Masa Örtüsü	150X200
10	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


HIP SURGERY PACK (WITH ELASTIC FENESTRATION)-1 /
KALÇA CERRAHİ SETİ (ELASTİK AÇIKLIK)-1



REF SV-1048

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	2	Side Banded Drape / Yan Bantlı Örtü	75x75
4	1	Hip Drape With Elastic Fenestration / Kalça Örtüsü Elastik Açıklık	240x300
5	1	Leg Cover / Bacak Kılıfı	40x120
6	1	Waste Bag / Atık Poşeti	28x33
7	2	Op Tape / Yapışkanlı Bant	10x50
8	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
9	1	Back Table Cover / Alet Masa Örtüsü	150x200


HIP SURGERY PACK (ELASTIC FENESTRATION WITH POCKET) /
KALÇA CERRAHİ SETİ (ELASTİK AÇIKLIK CEPLİ)



REF SV-1049

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	2	Side Banded Drape / Yan Bantlı Örtü	75x75
4	1	Hip Drape With Elastic Fenestration Double Side With Pocket / Kalça Örtüsü Elastik Açıklık Çift Taraf Cepli	300/240x350
5	1	Leg Cover (Stockinette) / Bacak Kılıfı	40x120
6	1	Waste Bag / Atık Poşeti	28x33
7	2	Op Tape / Yapışkanlı Bant	10x50
8	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
9	1	Back Table Cover / Alet Masa Örtüsü	150x200


SHOULDER ARTHROSCOPY PACK WITH POUCH /
OMUZ ARTROSKOPİ SETİ POŞLU



REF SV-1050

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	2	Banded Drape / Yan Bantlı Örtü	100x120
4	1	"U" Split Drape PE / "U" Yarıklı Örtü PE	150X200
5	1	"U" Split Drape With Pouch / "U" Yarıklı Poşlu omuz Örtüsü	200X300
6	1	Foot Drape / Ayak Örtüsü Bantlı	150x180
7	1	Stockinette / Stokinet	30x60
8	1	Elastic Bandage / Elastik Bandaj	ST
9	1	Skin Marker And Ruler / Cetvel ve Cilt Kalem	ST
10	4	Op-Tape / Yapıştırma Bandı	10x30
11	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
12	1	Back Table Cover / Alet Masa Örtüsü	150x200
13	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100


UPPER EXTREMITY PACK / ÜST EXTREMİTE SETİ



REF SV-1051

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	2	X-Ray Casette Bag / Röntgen Kaset Torbası	60X80
4	1	Foot Drape / Ayak Örtüsü Bantlı	150x180
5	1	Leg Cover Stockinette / Stokinet Bacak Kılıfı	30x60
6	1	Arm Drape / El Altı Örtüsü(Düz)	100x150
7	10	Gauze Swab / Gaz Spanç	10x10
8	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
9	1	Elastic Bandage / Elastik Bandaj	ST
10	1	Op-Tape / Yapıştırma Bandı	10x50
11	1	Esmach Bandage / Esmach Bandaj	5x200
12	1	Upper Extremity Drape / Üst Extremitte Örtüsü	200x240
13	1	Skin Marker And Ruler / Cetvel ve Cilt Kalemi	ST
14	1	Back Table Cover / Alet Masa Örtüsü	150x200
15	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100

VERTICAL ISOLATION DRAPE PACK /
VERTİKAL İZOLASYON ÖRTÜ SETİ



REF SV-1052

NO	QTY	NAME	SIZE
1	4	Towel / Havlu	40x40
2	2	Op Tape / Yapışkanlı Bant	10x50
3	1	Back Table Cover / Alet Masa Örtüsü	150X200
4	1	Vertical Isolation Drape Transparent PE (30*75 Incision Area) / Vertikal İzolasyon Örtüsü Şeffaf PE (30*75 iobanlı insz.)	240x320
5	1	Mayo Stand Cover / Mayo Masa Örtüsü	78x145
6	1	Waste Bag / Atık Poşeti	28x33


HAND SURGERY PACK / EL CERRAHİ ÖRTÜ SETİ



REF SV-1053

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	1	Hand Surgery Drape with Elastic Fenestration T Shaped / El Cerrahi Örtüsü (Elastik açıklık JT)	200/150x370
4	1	Stockinette / Stokinet	40x75
5	1	Elastic Bandage / Elastik Bandaj	10x150
6	1	Op Tape / Yapışkan Bant	10x50
7	1	Back Table Cover / Alet Masa Örtüsü	150x200
8	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100

PERCUTANEOUS PACK / PERKÜTAN SETİ



REF SV-1054

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	1	Percutaneous Drape / Perkütan Örtüsü	150x300
4	1	Scopy Cover / Skopi kılıfı	50x80
5	1	Camera Cover / Kamera Kılıfı	14x250
6	1	Op Tape / Yapışkanlı Bant	10x50
7	1	Back Table Cover / Alet Masa Örtüsü	150x200
8	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100

CYSTOSCOPY PACK / SİSTOSKOPİ SETİ



REF SV-1055

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	2	Leg Cover [Stockinette] / Bacak Kılıfı	60x130
4	1	Cystoscopy Drape / Sistoskopi Örtüsü	75x150
5	1	Basic Drape / Düz Örtü	75x90
6	1	Op Tape / Yapışkanlı Bant	10x30
7	1	Back Table Cover / Alet Masa Örtüsü	150x200
8	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100

TUR PACK / TUR SETİ



REF SV-1056

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	1	TUR Drape With Pouch / TUR Örtüsü L.T.T	175x240
4	1	Basic Drape / Düz Örtü	150x200
5	1	Back Table Cover / Alet Masa Örtüsü	100x150
6	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100

UROLOGY PACK / ÜROLOJİ SETİ



REF SV-1057

NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
3	1	Basic Drape / Düz Örtü	160x200
4	1	Urology Drape / Üroloji Örtüsü	175x240
5	1	Back Table Cover / Alet Masa Örtüsü	100x150
6	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100

DRESSING PACK / CERRAHİ PANSUMAN SETİ



REF SV-1058

NO	QTY	NAME	SIZE
1	1	Towel / Havlu	50x50
2	1	Gloves (Pair) / Eldiven (Çift)	7 Nr.
3	1	Tray / Tepsi	21x14x4.5
4	4	Gauze Swab / Gaz Spanç	7.5x7.5
5	2	Gauze Ball / Gaz Top	ST
6	1	Forceps / Forseps	ST
7	1	Cautery Forceps / Koter Forseps	ST


SUTURE PACK / SÜTUR SETİ



REF SV-1059

NO	QTY	NAME	SIZE
1	1	Basic Drape / Düz Örtü	50x50
2	1	Gloves (Pair) / Eldiven (Çift)	7-7.5 Nr.
3	1	Tray / Tepsi	21x14x4.5
4	4	Gauze Swab / Gaz Spanç	7.5x7.5
5	4	Gauze Ball / Gaz Top	ST
6	1	Surgical Blade / Bisturi	ST
7	1	Forceps / Forseps	ST
8	1	Cautery Forceps / Koter Forseps	ST
9	1	Waste Bag / Atık Torbası	28x33


DENTAL MINI PACK / DENTAL MİNİSETİ



REF SV-1060

NO	QTY	NAME	SIZE
1	1	Patient Gown / Bantlı Hasta Önlüğü	65x65
2	1	Transparent Tube Cover / Şeffaf Hortum Kılıfı	10x200
3	1	Transparent Tube Cover / Şeffaf Hortum Kılıfı	10x120
4	1	Reflector Drape / Reflektör Örtüsü	15x15
5	1	Table Drape / Tabla Örtüsü	40x45
6	1	Sterilization Wrap / Sterilizasyon Bohçası	50x50


IMPLANT PACK / İMPLANT SETİ



REF SV-1061


NO	QTY	NAME	SIZE
1	2	Standard Surgical Gown / Standart Cerrahi Önlük	L
2	1	Patient Gown / Bantlı Hasta Önlüğü	65x110
3	1	Transparent Tube Cover / Şeffaf Hortum Kılıfı	10x200
4	10	Gauze Swab / Gaz Spanç	7,5X7,5
5	2	Reflector Drape / Reflektör Örtüsü	15x15
6	1	Table Drape / Tabla Örtüsü	50x50
7	3	Bonnet / Bone	ST
8	2	Back Table Cover / Alet Masa Örtüsü	50x65
9	1	Sterilization Wrap / Sterilizasyon Bohçası	75x75

CHEMOTHERAPY PACK / KEMOTERAPİ SETİ



REF SV-1062


NO	QTY	NAME	SIZE
1	2	Towel / Havlu	40x40
2	1	Reinforced Chemotherapy Gown / Korumalı Kemoterapi Önlük	L
3	2	Chemotherapy Cloves / Kemoterapi Eldiveni	M
4	1	Chemotherapy Glasses / Kemoterapi Gözlüğü	ST
5	1	Face Mask / Maske	FFP3
6	1	Waste Bag / Tıbbi Atık Poşeti	50x80
7	1	Bonnet / Bone	ST
8	1	Plain Drape / Düz Örtü	40x80
9	1	Sterilization Wrap / Sterilizasyon Bohçası	75x 75

PLASTIC SURGERY PACK /
PLASTİK CERRAHİ SETİ

REF SV-1063

NO	QTY	NAME	SIZE
1	2	Side Banded Drape / Yan Bantlı Örtü	80x90
2	1	Foot Drape / Bantlı Ayak Örtüsü	150x180
3	1	Anesthesia Drape / Anestezi Örtüsü	150x240
4	1	Back Table Cover / Alet Masa Örtüsü	150x200
5	1	Sterilization Wrap / Sterilizasyon Bohçası	50x50


HAIR- IMPLANT PACK / SAÇ EKİM SETİ



REF SV-1064

NO	QTY	NAME	SIZE
1	1	Fenestrated Drape (10x19 With Adhesive) / Delikli Örtü (10x19 Bantlı)	75x75
2	2	Side Banded Drape / Yan Bantlı Örtü	75x75
3	1	Basic Drape / Düz Örtü	75x75
4	1	Sterilization Wrap / Sterilizasyon Bohçası	50x50


THORAX PACK / TORAKS SETİ



REF SV-1065

NO	QTY	NAME	SIZE
1	8	Towel / Havlu	40x40
2	2	Reinforced Surgical Gown / Korumalı cerrahi önlük	L
3	1	Side Banded Drape / Yan Bantlı Örtü	100x120
4	1	Thorax Drape With Incise Film (35X45) / Toraks Örtüsü (35X45 İnsizyonlu)	250/200x350
5	1	Op Tape / Yapıştırma Bandı	10x50
6	20	Gauze Swab / Gaz Spanç	7.5x7.5
7	1	Back Table Cover / Alet Masa Örtüsü	150x200
8	1	Sterilization Wrap / Sterilizasyon Bohçası	100x100



VETERINARY PACK / VETERİNER SETİ



REF SV-1066


NO	QTY	NAME	SIZE
1	1	Reinforced Surgical Gown / Korumalı Cerrahi Önlük	L
2	1	Basic Drape / Düz Örtü	120x150
3	1	Side Banded Drape / Yan Bantlı Örtü	100x150
4	1	Fenestrated Drape / Delikli Örtü	100x100
5	1	Back Table Cover / Alet Masa Örtüsü	100x120

CATHETERIZATION PACK / KATETER SETİ




REF SV-1067


NO	QTY	NAME	SIZE
1	1	Basic Drape / Düz Örtü	35x40
2	1	Fenestrated Drape [4 5x4 5 Fenestration] / Delikli Örtü [4 5X4 5 dık]	40x45
3	3	Gauze Swab / Gaz spanç	10x10
4	3	Gauze Ball / Gaz Top	ST
5	1	Syringe / Enjektör	10 cc
6	1	Forceps / Forseps	M
7	1	Cautery Forceps / Koter Forseps	ST
8	1	Gloves / Eldiven	ST
9	1	Gel / Jel	20 Gr.
10	1	Tray / Tepsi	21x14x4.5




Craniotomy Drape / **Kraniyotomi Örtüsü**
230x330




Carpal Tunnel Drape / **Karpal Tünel Örtüsü**
150x150




Spinal Vertebra Drape With Incise Film /
Spinal Vertebra Örtüsü (İnsizyon Film)
150x300




Shunt Drape / **Shunt Örtüsü**
150x300




Laparotomy Drape With Incise Film /
Laparotomi Örtüsü İnsizyonlu
240x200x320




Laparotomy Drape With Incise Film And Pouch /
Laparotomi Örtüsü İnsizyonlu-Poşlu
240x200x320




Cardiovascular Drape With "U" Split /
U- Yarıklı Kvc Örtüsü
240x200x340




Percutaneous Drape / **Perkütan Örtü**
150x300




Anesthesia Drape / **Anestezi Örtüsü**
150x240




Angiography Drape / **Anjiyo Örtüsü**
150x300




Angiography Drape With /
Panelli Anjiyo Örtüsü
195x300




Femoral Radial Angiography Drape With Double
Transparent Panels / **Femoral- Radyal Anjiyografi**
Örtüsü Çift Panelli
240x300




Arthroscopy Drape With Pouch /
Artroskopji Örtüsü Likit Torbalı
200x300




Hip Drape With Elastic Fenestration /
Kalça Örtüsü Elastik Film
240x300



Vertical Isolation Drape Transparent +Pe /
Vertikal İzolasyon Örtüsü Şeffaf +Pe
240x320




Ophthalmic Drape With Double Pouch /
Çift Poşlu Oftalmik Örtü
150x200




REF SV-1084

U Split Perineum Drape /
U - Yarıklı Perine Örtüsü
150x260




REF SV-1085

Hip Drape /
Kalça Örtüsü Çift Taraf Cepli
300x240x350




REF SV-1086

Urology Drape / Üroloji Örtüsü
150x300




REF SV-1087

Hand Surgery Drape
With Elastic Fenestration /
El Cerrahi Örtüsü
200x150x370




REF SV-1088

Arthroscopy Drape /
Artroskopi Örtüsü
200x300




REF SV-1089

Upper-Extremity Drape /
Extremite Örtüsü-Üst
150x240




REF SV-1090

Head Drape (Turban) / Baş Örtüsü
100x120



REF SV-1091

Gynaecology Drape With Pouch /
Jinekoloji Örtüsü-Poşlu
100x150




REF SV-1092

Gynaecology Drape / Jinekoloji Örtüsü
75x120




REF SV-1093

Under Buttocks Drape With Pouch /
Hasta Alt Örtüsü Poşlu
100x130




REF SV-1094

E.N.T Drape / KBB Örtüsü
160x240




REF SV-1095

O.P.U Drape Fenestrated /
O.P.U. Örtüsü
135x150



REF SV-1096

Fenestrated Drape / Delikli Örtü
80x100



REF SV-1097

Drape With Adhesive Fenestration /
Delikli Bantlı Örtü
100x100




REF SV-1098

Back Table Cover / Alet Masa Örtüsü
150x200




REF SV-1099

Basic Surgery Drape / Ameliyat Örtüsü
150x200




REF SV-1100

Food Drape / Ayak Örtüsü
160x180




REF SV-1101

Side Banded Drape / Yan Bantlı Örtü
100x120




REF SV-1102

Cardiovascular Drape / KVC Örtüsü
240x200x340




REF SV-1103

TUR Drape With Pouch /
Likit Toplamalı TUR Örtüsü
175x240




REF SV-1104

Laparoscopic Abdominal Perineal Drape /
Laparoskopik Abdominal Perinel Örtü
180x240x260




REF SV-1105

Laparoscopic Pelviscopy Drape /
Laparoskopi Pelviskopi Örtüsü
180x240x260



REF SV-1106

Thorax Drape With Incise Film /
Toraks Örtüsü
250x200x350




REF SV-1107

Thyroid Drape / Troid Örtüsü
240x320




REF SV-1108

Litotomic Gynaecology Drape With Pouch /
Litotomik Poşlu Jinekoloji Örtüsü
180x240x240




REF SV-1109

Litotomic Gynaecology Drape /
Litotomik Jinekoloji Örtüsü
180x240x240




REF SV-1110

Ophthalmic Drape With Single Pouch /
Tek Poşlu Oftalmik Örtü
100x130




REF SV-1111

Sub-Extremity Drape / Extremite Örtüsü-Alt
200x310



REF SV-1112

U-Split Drape +Pe / U-Yarıklı Örtü +Pe
240x350




REF SV-1113

U-Split Shoulder Arthroscopy Drape With Pouch /
U-Yarıklı Poşlu Omuz Artroskopik Örtüsü
240x320




REF SV-1114

Basic Surgery Drapes / Düz Örtü
150x250




REF SV-1115

Cystoscopy Drape / Sistoskopik Örtüsü
75x150



REF SV-1116

Mayo Stand Cover / Mayo Masa Örtüsü
90x140



REF SV-1117

Caesarean Drape / Sezeryan Örtüsü
200x300



DISPOSABLE VISITOR GOWN

Visitor gowns are used for minimal risk situations which contain slight barrier to small fluid penetration. The gowns are breathable and comfortable to wear.



SIZE	XS	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	56	56	56	56	56	56	58	58	+/- 3 cm
B	58	59	60	61	62	63	64	65	+/- 3 cm
C	105	110	115	120	125	130	135	140	+/- 3 cm
D	130	135	140	145	150	155	160	165	+/- 3 cm

SV- 180

DETAILS

Model: SV-180 series
Fabric: PP (spunbond+spunbond)
Weight: 17 - 80 GSM
Cuff: Elastic/Knitted
Color: White/Yellow/Blue/Green
Stitch: Overlock
Other: Waist and neck straps for a secure fit
Optional: Tie on neck - Velcro neck
Sterile/Non Sterile

For more informatin
please refer to the data sheets.



DISPOSABLE ISOLATION GOWN (PP+PE)

DISPOSABLE SURGICAL GOWN (SMMS-SSMMS)

HIGH PROTECTION SURGICAL GOWN (SSMLMS + PE)

SV- 140

DETAILS

Model: SV- 140 series
Fabric: PP + PE (polyethylene) laminated
Weight: 30-80 GSM
Cuff: Elastic/Knitted
Color: Blue/Yellow/Green/White
Stitch: Overlock / Ultrasonic seam
Other: Waist and neck straps for a secure fit
Tie on neck - Velcro neck
Sterile/Non Sterile

For more informatin
please refer to the data sheets.

SV- 141-142

DETAILS

Model: SV- 141 - 142 series
Fabric: SMMS - SSMMS
Weight: 25-80 GSM
Cuff: Elastic/Knitted
Color: White/Yellow/Blue/Green
Stitch: Overlock/Ultrasonic
Other: Waist and neck straps for a secure fit
Optional: Sterile/Non Sterile Velcro neck

For more informatin
please refer to the data sheets.

SV- 143

DETAILS

Model: SV- 143 series
Fabric: 17 gsm SSMMS Nonwoven + 1.5 gsm Hotmelt + 18 gsm Triple-layered Breathable Medical PE Film + 1.5 gsm Hotmelt + 17 gsm SSMMS Nonwoven
Weight: 55-90 GSM
Cuff: Elastic/Knitted
Color: Blue
Stitch: Overlock/ Ultrasonic seam
Other: Waist and neck straps for a secure fit
Optional: Static/Antistatic
Tie on neck - Velcro neck
Sterile/Non Sterile

For more informatin
please refer to the data sheets.



SIZE	XS	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	56	56	56	56	56	56	58	58	+/- 3 cm
B	58	59	60	61	62	63	64	65	+/- 3 cm
C	105	110	115	120	125	130	135	140	+/- 3 cm
D	130	135	140	145	150	155	160	165	+/- 3 cm

COLORS



SERVICES





DISPOSABLE **Reinforced**
SURGICAL GOWN (SMS+PE)

SV- 144

DETAILS

Model: SV- 144 series
Fabric: SMS + PE (spunbond+meltblown
+spunbond + polyethylene
laminated)
Weight: 45-80 GSM
Cuff: Elastic/Knitted
Color: Blue
Stitch: Overlock/ Ultrasonic seam
Other: Waist and neck straps for a secure fit
Optional: Tie on neck - Velcro neck
Sterile/Non Sterile

For more informatin
please refer to the data sheets.

DISPOSABLE **Reinforced**
SURGICAL GOWN / **Ultimate**

SV- 144

DETAILS

Model: SV- 144 series
Fabric: SMS + PE (spunbond+meltblown
+spunbond + polyethylene
laminated)
Weight: 45-80 GSM
Cuff: Elastic/Knitted
Color: Blue
Stitch: Overlock/ Ultrasonic seam
Other: Waist and neck straps for a secure fit
Optional: Tie on neck - Velcro neck
Sterile/Non Sterile

For more informatin
please refer to the data sheets.



SIZE

	XS	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	56	56	56	56	56	56	58	58	+/- 3 cm
B	58	59	60	61	62	63	64	65	+/- 3 cm
C	105	110	115	120	125	130	135	140	+/- 3 cm
D	130	135	140	145	150	155	160	165	+/- 3 cm

COLORS



SERVICES






DISPOSABLE EXAMINATION GOWN (SMS)

SV- 145

DETAILS
Model: SV- 145 series
Fabric: SMS (spunbond+meltblown+spunbond)
Weight: 20-80 GSM
Cuff: Elastic/Knitted
Color: White/Yellow/Blue/Green
Stitch: Overlock - Ultrasonic
Other: Waist and neck straps for a secure fit
Optional: Tie on neck - Velcro neck
Sterile/Non Sterile

For more informatin
please refer to the data sheets.




DISPOSABLE SURGICAL GOWN (Thermoband Striped Gown)

SV- 146

DETAILS
Model: SV- 146 series
Fabric: PP + PE + Thermoband
Weight: 40 - 80 GSM
Cuff: Elastic/Knitted
Color: Blue/Yellow/Green/White
Stitch: Overlock/Welded
Type: PB-3B / PB-4B / PB- 6B AAMI Level 2 - 3
Other: Waist and neck straps for a secure fit
Optional: Sterile/Non Sterile
Velcro neck
Ultrasonic seam

For more informatin
please refer to the data sheets.



DISPOSABLE SURGICAL GOWN / Classic

SV- 147

DETAILS
Model: SV- 147 series
Fabric: Spunlace / Sontara
Weight: 40 - 95 GSM
Cuff: Elastic/Knitted
Color: Blue/White
Stitch: Overlock/Welded
Type: PB-3B / PB-4B / PB- 6B AAMI Level 2 - 3
Other: Waist and neck straps for a secure fit
Optional: Sterile/Non Sterile
Velcro neck
Ultrasonic seam

For more informatin
please refer to the data sheets.



SIZE	XS	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	56	56	56	56	56	56	58	58	+/- 3 cm
B	58	59	60	61	62	63	64	65	+/- 3 cm
C	105	110	115	120	125	130	135	140	+/- 3 cm
D	130	135	140	145	150	155	160	165	+/- 3 cm

SERVICES





V-NECK SET

SV- 120

DETAILS

Model: SV- 120 series
Fabric: SS-SMS-SMMS
Weight: 30-80 GSM
Color: Blue/Yellow/Green/White
Stitch: Overlock/Welded
Other: Desings with two or three pockets
Optional: Sterile/Non Sterile

For more informatin
please refer to the data sheets.

THERMAL SET

SV- 121

DETAILS

Model: SV- 121 series
Fabric: SS-SMS-SMMS
Weight: 30-80 GSM
Cuff: Knitted
Color: Blue/Yellow/Green/White
Stitch: Overlock/Welded
Other: Desings with two or three pockets
Optional: Sterile/Non Sterile

For more informatin
please refer to the data sheets.

THERMAL JACKET

SV- 122


DETAILS

Model: SV- 122 series
Fabric: SS-SMS-SMMS
Weight: 30-80 GSM
Cuff: Knitted
Color: Blue/Yellow/Green/White
Stitch: Overlock/Welded
Other: Desings with two or three pockets
Optional: Sterile/Non Sterile

For more informatin
please refer to the data sheets.

		S	M	L	XL	2XL	3XL	Average
	A	74	75	76	77	78	79	+/- 3cm
	B	60	62	64	66	68	70	+/- 3cm
	C	35	36	37	38	39	40	+/- 3cm
	D	18	19	20	21	22	23	+/- 3cm

		S	M	L	XL	2XL	3XL	Average
	A	102	104	106	108	110	112	+/- 3cm
	B	52	54	57	61	64	67	+/- 3cm
	C	54	56	59	63	66	69	+/- 3cm
	D	21	22	23	24	25	26	+/- 3cm



		S	M	L	XL	2XL	3XL	Average
	A	73	75	77	79	81	83	+/- 3cm
	B	55	57	59	61	63	65	+/- 3cm
	C	59	60	61	62	63	64	+/- 3cm
	D	28	28	28	28	29	29	+/- 3cm

		S	M	L	XL	2XL	3XL	Average
	A	102	104	106	108	110	112	+/- 3cm
	B	52	54	57	61	64	67	+/- 3cm
	C	56	56	59	63	66	69	+/- 3cm
	D	21	22	23	24	25	26	+/- 3cm



KIMONO SET

SV- 123

DETAILS

Model: SV- 123 series
Fabric: SMS (spunbond+meltblown+spunbond)
Weight: 20-80 GSM
Color: White/Yellow/Blue/Green
Stitch: Overlock - Ultrasonic
Other: Waist and neck straps for a secure fit
Optional: Tie on neck - Velcro neck
Sterile/Non Sterile

For more information
please refer to the data sheets.

COLORS



SERVICES




LAB COAT

SV- 189

DETAILS

Model: SV- 189 series
Fabric: SS-SMS-SMMS
Weight: 20-80 GSM
Cuff: Elastic, Knitted
Color: Blue/Yellow/Green/White
Stitch: Overlock/Welded
Other: Designs with two or three pockets
Optional: Sterile/Non Sterile

For more information
please refer to the data sheets.



	S	M	L	XL	2XL	3XL	4XL	Average
A	78	81	84	87	90	93	96	+/- 3cm
B	125	130	135	140	145	150	155	+/- 3cm
C	26	27	28	29	30	31	32	+/- 3cm

	S	M	L	XL	2XL	3XL	4XL	Average
A	102	104	106	108	110	112	114	+/- 3cm
B	52	54	57	61	64	67	70	+/- 3cm
C	54	56	59	63	66	69	72	+/- 3cm
D	21	22	23	24	25	26	27	+/- 3cm



	S	M	L	XL	2XL	3XL	4XL	Average
A	96	99	102	105	108	111	114	+/- 3cm
B	61	64	67	70	73	76	79	+/- 3cm
C	56	56	56	56	56	58	58	+/- 3cm



DISPOSABLE COVERALL TYPE 5B/6B

DISPOSABLE COVERALL TYPE 3B/4B THERMO BAND STRIPED

DISPOSABLE COVERALL TYPE 3B/4B ULTRASONIC

SV- 100

DETAILS

Model: SV- 100 series

Fabric: Spunbond hotmelt polyethylen, PP + PE

Weight: 40-80 GSM

Cuff: Elastic

Color: White

Stitch: Overlock

Other: Waist and neck straps for a secure fit

Optional: Static/Antistatic

Optional: Sterile/Non Sterile

Optional: Ultrasonic seam

For more informatin please refer to the data sheets.

SV- 110

DETAILS

Model: SV- 110 series

Fabric: Spunbond hotmelt polyethylen, PP + PE

Weight: 55-80 GSM

Cuff: Elastic

Color: White

Stitch: Overlock

Other: Waist and neck straps for a secure fit

Optional: Static/Antistatic

Optional: Sterile/Non Sterile

For more informatin please refer to the data sheets.

SV- 113

DETAILS

Model: SV- 113 series

Fabric: Spunbond hotmelt polyethylen, PP + PE

Weight: 40-80 GSM

Cuff: Elastic

Color: White

Stitch: Ultrasonic

Other: Waist and neck straps for a secure fit

Optional: Static/Antistatic

Optional: Sterile/Non Sterile

For more informatin please refer to the data sheets.

SIZE	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	55	60	65	70	75	80	85	+/- 3 cm
B	188	190	192	194	196	198	200	+/- 3 cm
C	74	76	78	80	82	84	86	+/- 3 cm

COLORS



COLORS



SERVICES





DISPOSABLE COVERALL
(SMS - PP+PE)
Without Hood


DISPOSABLE COVERALL
(SMS)

HI-VIS REFLECTIVE
WORKWEAR COVERALL

SV- 107

DETAILS
Model: SV- 107 series
Fabric: SMS, PP+PE
Weight: 20-80 GSM
Cuff: Elastic / Knitted cuff
Color: Blue,Yellow,Green,White
Stitch: Overlock
Optional: Static/Antistatic


For more informatin
please refer to the data sheets.



SV- 188

DETAILS
Model: SV-188 series
Fabric: Spunbond+meltblown+ spunbond
Weight: 20-80 GSM
Cuff: Elastic / Knitted cuff
Color: Blue,Yellow,Green,White
Stitch: Overlock


For more informatin
please refer to the data sheets.




SV- 116

DETAILS
Model: SV- 116 series
Fabric: SMS or Microporous PE Laminated
Weight: 40-90 GSM
Cuff: Elastic
Stitch: Overlock/ Ultrasonic stripe
Ergonomically shaped hood
Reflective tape
Anti-static
Optional: Sterile/Non Sterile

For more informatin
please refer to the data sheets.



SIZE	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	55	60	65	70	75	80	85	+/- 3 cm
B	160	162	164	166	168	170	172	+/- 3 cm
C	74	76	78	80	82	84	86	+/- 3 cm



SIZE	S	M	L	XL	2XL	3XL	4XL	Tolerans
A	55	60	65	70	75	80	85	+/- 3 cm
B	188	190	192	194	196	198	200	+/- 3 cm
C	74	76	78	80	82	84	86	+/- 3 cm



FFP2 MASKS
FFP2 MASKS NR D

SV- 130 SV- 131
SV- 130-D SV- 131-D

DETAILS

Model: SV- 130-SV-131
SV-130-D - SV-131-D
Fabric: Spunbond + Meltblown
Sizes: Small / Medium / Standard
Weight: 150-180 GSM
Color: Black/White
Filtration: Efficiency ≥ 95%

For more informatin
please refer to the data sheets.

FFP3 MASKS
FFP3 MASKS NR D

SV-200 SV-200D

DETAILS

Model: SV-200 - SV-200D
Fabric: Spunbond + Meltblown
Sizes: Small / Medium / Standard
Weight: 150-180 GSM
Color: White
Filtration: Efficiency ≥ 95%

For more informatin
please refer to the data sheets.

SURGICAL MASKS


2000R

DETAILS

Fabric: Spunbond + Meltblown + Spunbond
Weight: 70-80 GSM
Stitch: Ultrasonic seam
Filtration: BFE > 95%
Other: EU Type IIR I ASTM LVL 1

Adjustable nose wire , Anti Allergic
Air Permeable, Soft ear loops
Splash resistant, Hygienic and hypoallergenic

For more informatin
please refer to the data sheets.





SV- 49 | Round Elastic Bonnet Cap

DETAILS

Fabric : SS , SMS, PP+PE Double stitch
Weight : 20- 60 gr Anti_allergic.



SV- 50 | Operation Bonnet

DETAILS

Fabric : SS , SMS Double stitch
Weight : 20- 60 gr Anti_allergic.



SV- 51 | Astro Hood

DETAILS

Fabric : SS , SMS Double stitch
Weight : 20- 60 gr Anti_allergic.



SV- 52 | Barret Hood

DETAILS

Fabric : SS Elastic
Weight : 10- 20 gr Anti_allergic.



SV- 53 | Surgical Hoods With Ties

DETAILS

Fabric : SS , SMS, PP+PE Elastic
Weight : 20- 60 gr Anti_allergic.



SV- 54 | Headgear With Cape

DETAILS

Fabric : SS , SMS, PP+PE Elastic
Weight : 20- 60 gr Anti_allergic.



SV- 55 | PE Cap

DETAILS

Fabric : PE Elastic
Weight : 9 - 20 Micron Anti_allergic.



SV- 150 | Boot Cover

DETAILS

Fabric: PP +PE
Weight: 30 -70 GSM
Stitch: Overlock
Optional: Ultrasonic seams, Non-slip velcro tape



SV- 151 | Thermo Band Striped Long Boot Cover

DETAILS

Fabric: PP+PE + Thermoband
Weight: 40-80 GSM
Stitch: Overlock
Optional: Ultrasonic seams, Non-slip velcro tape



SV- 152 | Short Overshoe

DETAILS

Fabric: PP+PE
Weight: 30-80 GSM
Stitch: Overlock
Optional: Ultrasonic seams, Non-slip velcro tape



SV- 60 | Shoe Cover Simple

DETAILS

Fabric: PE Thickness: 9-40 Micron
Stitch: Rope stitch technology is used to avoid shoe cover ladders
Other: Reinforced double elastic band, liquid-proof, non-slip, tear-proof - Produced in automatic machines without being touched Compressed vacuum packaging




SV- 63 | SMS Long Boot Cover

DETAILS

Fabric: SS or SMS
Weight: 20-80 GSM
Stitch: Overlock
Optional: Ultrasonic seams, Non-slip velcro tape





Biodegradable



SV- 630 Hairdresser Towel	FOR WOMEN	SIZE 40 x 80 cm (16" x 31")	THICKNESS STANDART ☒ THICK ☐ 50gsm / m ²	1 PER PACK 1 DO NOT REUSE
SV- 631 Hairdresser Towel	FOR MEN	SIZE 30x60 cm (12" x 24")	THICKNESS STANDART ☒ THICK ☐ 50gsm / m ²	1 PER PACK 1 DO NOT REUSE
SV- 632 Feet Towel		SIZE 40 x 40 cm (16" x 16")	THICKNESS STANDART ☐ THICK ☒ 70gsm / m ²	1 PER PACK 1 DO NOT REUSE

SV- 633 | Hand Towel

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for travel, camping, long rides, or daily use.

SIZE
40 x 40 cm
(16" x 16")

1
PER
PACK

DO NOT
REUSE

THICKNESS
STANDART THICK
50gsm / m²

SV- 634 | Sports/Gym Towel

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for sport, travel, camping, long rides, or daily use.

SIZE
40 x 80 cm
(16" x 31")

1
PER
PACK

DO NOT
REUSE

THICKNESS
STANDART THICK
50gsm / m²

SV- 635 | Bath Towel

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for travel, camping, long rides, or daily use.

SIZE
60 x 120 cm
(24" x 47")

1
PER
PACK

DO NOT
REUSE

THICKNESS
STANDART THICK
50gsm / m²

SV- 636 | Bath Towel

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for sport, travel, camping, long rides, or daily use.

SIZE
80 x 140 cm
(31" x 55")

1
PER
PACK

DO NOT
REUSE

THICKNESS
STANDART THICK
50gsm / m²

SV- 637 | Pet Towel

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for Pet, travel, camping, long rides, or daily use.

SIZE
80 x 140 cm
(31" x 55")

1
PER
PACK

DO NOT
REUSE

THICKNESS
STANDART THICK
70gsm / m²



Disposable Towel Set

SV- 638 | Towel Set
Set of 3 (Hair, bath, foot)

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for travel, camping, long rides, or daily use.

SIZE / HAIR
40 x 80 cm
(16" x 31")

SIZE / BATH
80 x 140 cm
(31" x 55")

SIZE / FOOT
40 x 45 cm
(16" x 18")

THICKNESS
STANDART THICK
50gsm / m²

1
PER
PACK

DO NOT
REUSE



SV- 639 | Towel Set
Set of 2 (Hair, bath)

Absorbs water and dries body after showering
Resists tearing when wet for more drying Great for travel, camping, long rides, or daily use.

SIZE / HAIR
40 x 80 cm
(16" x 31")

SIZE / BATH
80 x 140 cm
(31" x 55")

THICKNESS
STANDART THICK
70gsm / m²

1
PER
PACK

DO NOT
REUSE

Disposable Bed Set

SV- 640 | Single Bed Set
1 Bed sheet, 1 pillowcase, 1 duvet cover

Disposable duvet cover set is hygienic and ultra-light. Travel, Hospital, hotel, camping, caravan, yacht, Personal use...

BED SHEET
160x 200 cm
(64" x 80")

PILLOW CASE
50x 70 cm
(20" x 28")

DUVET COVER
160x 200 cm
(64" x 80")

THICKNESS
STANDART THICK
50gsm / m²

1
SET
PACK

DO NOT
REUSE

SV- 641 | Double Bed Set
1 Bed sheet, 2 pillowcase, 1 duvet cover

Disposable duvet cover set is hygienic and ultra-light. Travel, Hospital, hotel, camping, caravan, yacht, Personal use...

BED SHEET
200x 240 cm
(80" x 96")

PILLOW CASE
50x 70 cm
(20" x 28")

DUVET COVER
160x 200 cm
(80" x 96")

THICKNESS
STANDART THICK
50gsm / m²

1
SET
PACK

DO NOT
REUSE



Biodegradable

SV- 649 | Spa Massage Sheet

Disposable duvet cover set is hygienic and ultra-light. Hospital, hotel, yacht, Professional use...

SIZE
80 x 200 cm
(31" x 79")

1 PER PACK

DO NOT REUSE

THICKNESS

STANDARD THICK

50gsm / m²

SV- 650 | Spa Massage Sheet (Hole)

Disposable duvet cover set is hygienic and ultra-light. Hospital, hotel, yacht, Professional use...

SIZE
80 x 200 cm
(31" x 79")

1 PER PACK

DO NOT REUSE

THICKNESS

STANDARD THICK

50gsm / m

SV- 651 | Sun Lounger Cover

Disposable duvet cover set is hygienic and ultra-light. Hospital, hotel, yacht, Professional use...

SIZE
80 x 220 cm
(31" x 87")

1 PER PACK

DO NOT REUSE

THICKNESS

STANDARD THICK

50gsm / m²

SV- 652 | Sun Lounger Cover With Pouch

Disposable duvet cover set is hygienic and ultra-light. Hospital, hotel, yacht, Professional use...

SIZE
80 x 200 cm
(31" x 79")

1 PER PACK

DO NOT REUSE

THICKNESS

STANDARD THICK

50gsm / m

Disposable Travel Set

SV- 642 | Single Travel Set

1 Bath towel, 1 Hair towel, 1 Foot towel
1 Bed Sheet, 1 Pillow Case

Disposable duvet cover set is hygienic and ultra-light. Travel, Hospital, hotel, camping, caravan, yacht, Personal use...

THICKNESS

STANDARD THICK

50gsm / m²

1 SET PACK

DO NOT REUSE

BATH TOWEL
80x 140 cm
(31" x 55")

FOOT TOWEL
40x 45 cm
(16" x 18")

HAIR TOWEL
40x 80 cm
(16" x 31")

PILLOW CASE
50x 70 cm
(20" x 28")

DUVET COVER
160x 200 cm
(63" x 79")

SV- 643 | Double Travel Set

2 Bath towel, 2 Hair towel, 2 Foot towel
1 Bed Sheet, 2 Pillow Case

Disposable duvet cover set is hygienic and ultra-light. Travel, Hospital, hotel, camping, caravan, yacht, Personal use...

THICKNESS

STANDARD THICK

50gsm / m²

1 SET PACK

DO NOT REUSE

BATH TOWEL
80x 140 cm
(31" x 55")

FOOT TOWEL
40x 45 cm
(16" x 18")

HAIR TOWEL
40x 80 cm
(16" x 31")

PILLOW CASE
50x 70 cm
(20" x 28")

DUVET COVER
160x 200 cm
(63" x 79")



SV- 644 | Double Duvet Case

Disposable duvet cover set is hygienic and ultra-light. Travel, hospital, hotel, camping, caravan, yacht, Personal use...

SIZE
200 x 240 cm
(79" x 96")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 PER PACK 1 DO NOT REUSE

SV- 645 | Single Duvet Case

Disposable duvet cover set is hygienic and ultra-light. Travel, hospital, hotel, camping, caravan, yacht, Personal use...

SIZE
160 x 200 cm
(63" x 79")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 PER PACK 1 DO NOT REUSE

SV- 646 | Double Bed Sheet

Disposable duvet cover set is hygienic and ultra-light. Travel, hospital, hotel, camping, caravan, yacht, Personal use...

SIZE
200 x 240 cm
(79" x 96")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 PER PACK 1 DO NOT REUSE

SV- 647 | Single Bed Sheet

Disposable duvet cover set is hygienic and ultra-light. Travel, hospital, hotel, camping, caravan, yacht, Personal use...

SIZE
160 x 200 cm
(63" x 79")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 PER PACK 1 DO NOT REUSE

SV- 648 | Single Pillow Case

Disposable duvet cover set is hygienic and ultra-light. Travel, hospital, hotel, camping, caravan, yacht, Personal use...

SIZE
50 x 80 cm
(20" x 31")

THICKNESS
STANDART ☒ THICK ☐
50gsm / m²

1 PER PACK 1 DO NOT REUSE

Disposable Bed Sheet



BED SHEET

Bed covers are ideal for any type of medical table, stretcher or patient bed. This includes hospital bed mattress cover for ambulance stretchers through to hospitality linens. They can be used in patient rooms in hospitals, clinics, or patient care facilities.

DETAILS

Fabric: SMS, SS
Weight: 35-60 GSM
Color: White/Green/Blue

COLORS



BED COVER

Made of non-woven fabric, disposable bed cover can reduce the risk of cross-infection effectively. Bed cover is soft, economic, lightweight, and breathable. Bed covers are also equipped with elastic bands on the two sides which is more convenient to use, and offer more protection. It is one of the best choices for patients' health in medical facilities.

DETAILS

Fabric: SMS, SS
Weight: 35-60 GSM
Color: White/Green/Blue

SERVICES



EU-MDR

Certificates Modul B



Certificates Modul C2



Certificates




Test Reports





■ Raw Material



■ Tailoring



■ Sewing



■ Pressure Plastic Molding



■ Packaging



■ Delivery

How is production carried out

Production process and management are conducted based on the internationally accepted Strict hygiene standards. The products are duly packaged and delivered to the end users on time following the guiding hygiene principles.


 Fast delivery

Fast delivery from our Europe central warehouse with high capacity Turkey based production.

 EU certification

In order to be able to offer the best protection and the greatest security in everyday working life, we only supply certified products.

 Top products

We only carry tested, certified and high quality products.

 Fair prices

Thanks to strong partnerships and optimized business processes, we can deliver on our products at fair prices.

Declaration On Product Standards & Services



CATEGORY 3

Liquid-tight protective clothing: Full body or full protection suits against liquid chemicals



EN 13688

General requirement for protective clothing.
EN ISO 13688:2013



EN 14126

Protective clothing against infectious agents.
EN 14126: 2003 + AC:2004



EN 14605 TYPE PB 3-B

Liquid-tight protective clothing against liquid chemicals.
EN 14605:2005 + A1:2009 TYP PB 3-B



EN 13034 TYPE PB 6-B

Protective clothing against liquid chemicals with limited protection.
EN 13034:2005 + A1:2009 TYP PB 6-B



EN 14605 TYPE PB 4-B

Spray-tight protective clothing against liquid chemicals.
EN 14605:2005+A1:2009



EN 13982 TYPE 5-B

Protective clothing against airborne solid particles
EN ISO 13982-2:2004 TYP 5



EN 1149

Protective clothing with electrostatic properties
EN 1149-1:2006

EN 13795

Surgical drapes, gowns and clean air suits, used as medical devices for patients, clinical staff and equipment.



EN 1073-2:2002

Protective clothing against radioactive contamination-part2

FFP2 NR D

Passes EN 149:2001 optional Dolomite (D) clogging test for more comfortable breathing.

FFP3 NR D

Passes EN 149:2001 optional Dolomite (D) clogging test for more comfortable breathing.

Services



HOSPITAL



INDUSTRY



CORONA PROTECTION



EMERGENCY



HYGIENE



LABORATORY

PODİMA

SUVICOM®

rowimed®

Reddax®



About Podima

Our products are produced in accordance with our production infrastructure, expert staff and strict hygiene standards in Türkiye. Accordingly, our products are used and demanded safely in Europe and all over the world. Our quality, production and distribution performance has been appreciated by all our customers, especially during the pandemic period.

Our **SUVICOM** brand is primarily used in Germany and our other brand, **ROWIMED**, is used in Italy with its high quality and with confidence

Our Philosophy Our goal is to supply quality medical PPEs and supplies reliably and continuously at fair prices. We are committed to ensure a gap-free supply chain system, and achieve highest customer satisfaction and sustainable partnerships.

As a manufacturer aiming for the highest production standards, quality is an essential element for us and therefore we implement a rigorous control mechanism. Quality controls, certification and compliance, sustainable supply chain and customer satisfaction are our core principles. **As Podima Medikal**, we make constant effort to maintain these standards and prioritize quality at every stage of our supply chain.

Imprint / Contact

Türkiye
Güçlü
Pantolonlu
Pantolonlu

Head Office
Podima Medikal ve
Tekstil San. Tic. Ltd. Şti.
Şerifali Mh, Mevdudi Sk.
No:53/1, 34775 Ümraniye
İstanbul / Türkiye

Phone: +90 216 784 07 93
+90 216 314 93 16
Email: info@podimamedikal.com

■ Factory:
Siirt Yolu Cad. 72100
Merkez Batman / Türkiye

■ Branch Office
Suivicom Media GmbH
Lilienthalstr 14
41515, Grevenbroich,
Germany

Phone: +49 221 292033 31
Email: info@suivicom.de
Web: www.suivicom-medical.de